

Using AutoEncoder 5.0

Table of Contents

PREFACE.....	7
1.0 Website Terms and Conditions	7
2.0 Intended Audience.....	10
3.0 Manual Conventions	10
4.0 Getting Started.....	10
LOG IN/LOG OFF	11
4.1 Log In.....	11
4.2 Change Password	12
4.2.1 Password Rules	13
4.2.2 Retrieve a Forgotten Password	14
4.2.3 Password Notifications	15
4.2.4 Log Off.....	15
USING AUTOENCODER	17
5.0 About AutoEncoder.....	17
6.0 AutoEncoder Workflow.....	17
7.0 AutoEncoder Interface	19
Main AutoEncoder Window	19
7.1 Home Page (Dashboard)	19
7.2 Status Bar	20
7.3 Menus	21
Dashboard	24
7.4 Use the Dashboard for Coding.....	24
7.4.1 View Coding Stages	25
7.4.2 View a Term's Audit History	27
7.4.3 Confirm/Review/Approve Terms	28
7.4.4 Re-autocode Terms.....	30
7.4.5 Code/Unicode Matched Terms.....	31
7.4.6 View Autoencoding Job Details.....	33
7.4.7 View/Unicode Terms in Transmission Queue	34
7.4.8 View Module Reports	35
Basic AutoEncoder Functions	36
7.5 Update Your Account Information	36
7.6 View System Versions	38
7.7 Return to Home Page.....	39

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

7.8 Remove Message	39
7.9 Change Role	40
Table Functions.....	40
7.10 Use Table Functions.....	40
7.11 Sort and Filter Tables	41
7.11.1 Sort a Column	41
7.11.2 Move a Column.....	42
7.11.3 Show/Hide a Column	42
7.11.4 Filter a Column on Criteria.....	43
7.11.5 Resize a Column	43
7.12 Export to Excel	44
7.13 Use Information Icon	45
7.14 Refresh Data	46
7.15 Reset Table to Default Layout.....	47
OBJECTS	49
8.0 Manage Objects	49
8.1 Edit an Object.....	49
8.2 Delete an Object	50
8.3 Edit Object Synonym Lists	51
8.4 View Object Reports	52
8.5 Edit Object Settings.....	53
Trials	55
8.6 Manage Trials.....	55
8.6.1 View Trials.....	55
8.6.2 Add a New Trial	57
8.7 Mapping Overview	59
8.7.1 Manage Mapping	63
8.7.2 View/Edit a Mapping Between Modules and Dictionary Terms	63
8.7.3 Create a Mapping Between Modules and Dictionary Terms	65
8.7.4 Upgrade Dictionary Version for a Module's Mapping.....	67
8.7.5 Delete a Mapping Between Modules and Dictionary Terms	68
Trial Groups	70
8.8 Manage Trial Groups.....	70
8.8.1 View Trial Groups.....	70
8.8.2 Add a New Trial Group.....	71
8.8.3 View Trials within a Trial Group	72

Organizations.....	72
8.9 Manage Organizations	72
8.9.1 View Organizations	73
8.9.2 Add a New Organization	74
8.9.3 View Trials within an Organization	75
USERS.....	77
9.0 Manage Users	77
9.1 View Users	77
9.2 Edit/Add a User	78
9.3 Reset a User's Password	80
9.4 Assign User to an Object	81
9.5 View User Reports	82
9.6 User Permissions	83
ROLES.....	85
10.0 Manage Roles	85
10.1 View Roles.....	86
10.2 Assign Role to an Object	87
DICTIONARIES	89
11.0 Manage Dictionaries	89
11.1 View Dictionaries	89
11.2 Add New Dictionary Version	91
11.3 Delete Dictionary Version	92
SYNONYM LISTS	93
12.0 Manage Synonym Lists	93
12.1 View Synonym Lists	93
12.2 Edit/Add Synonym Lists	95
12.3 Delete a Synonym List.....	96
12.4 Import/Export Synonym List Information	97
12.5 Manage Synonym Terms	98
12.5.1 Add Synonym Terms	98
12.5.2 Delete Synonym Terms	100
SYSTEM REPORTS.....	101
13.0 View System Reports	101
SYSTEM SETTINGS.....	103
14.0 View/Edit System Settings	103

Revision History

Revision	Description of Change	Author	Effective Date
1.0	New document	Lorrie Tate	23-APR-2018

PREFACE

1.0 Website Terms and Conditions

YOUR USE AND ENJOYMENT OF OUR WEBSITE IS SUBJECT TO VARIOUS TERMS AND CONDITIONS SET FORTH BELOW. PLEASE READ THE FOLLOWING TERMS AND CONDITIONS CAREFULLY. YOUR USE OF THIS WEBSITE CONSTITUTES YOUR AGREEMENT TO BE BOUND BY ALL OF SUCH TERMS AND CONDITIONS, WITHOUT MODIFICATION, AS SUCH TERMS AND CONDITIONS MAY BE MODIFIED BY OMNICOMM SYSTEMS, INC. FROM TIME-TO-TIME AND POSTED ON THIS WEBSITE.

Website Terms and Conditions

1. Copyrights and Other Intellectual Property. Except as otherwise noted in this Website, all information, documentation, software, applications and other content posted in this Website (or accessed through this Website), including, without limitation, all text, graphics, logos, icons and images and the compilation, collection, arrangement, overall appearance and method of assembly of same (collectively, the “**Information**”), is the sole and exclusive property of OmniComm Systems, Inc. (the “**Company**”), or its third party vendors and content suppliers. The Information is protected by U.S. and international copyright laws. You may not sell, republish, copy for other than your own individual use, modify, distribute, license, create derivative works from or otherwise distribute the Information, or any portion thereof, without the prior written consent of the Company in each instance, which consent shall be given or withheld in the Company’s sole and absolute discretion.
2. **NO WARRANTIES; LIMITATION OF LIABILITY.** ALTHOUGH THE COMPANY ATTEMPTS TO MAKE SURE THE INFORMATION ON THIS WEBSITE IS ACCURATE AND CURRENT, THE COMPANY MAKES NO WARRANTIES OR REPRESENTATIONS OF ANY NATURE OR KIND THAT THE INFORMATION IS ACCURATE OR CURRENT OR THAT THE INFORMATION IS SUITABLE FOR ANY PARTICULAR PURPOSE. ALL INFORMATION ON THIS WEBSITE IS PROVIDED TO YOU “AS IS” WITHOUT WARRANTY OF ANY KIND, INCLUDING EXPRESS OR IMPLIED WARRANTIES OF MERCHANTABILITY, FITNESS FOR A PARTICULAR PURPOSE, TITLE AND NON-INFRINGEMENT. IN NO EVENT SHALL THE COMPANY AND/OR ITS RESPECTIVE VENDORS OR SUPPLIERS OR OTHER PARTIES INVOLVED IN CREATING OR DELIVERING THIS WEBSITE OR THE INFORMATION ON THIS WEBSITE BE LIABLE FOR ANY DAMAGES WHATSOEVER ARISING OUT OF OR RELATED TO YOUR USE OF OR ACCESS TO THIS WEBSITE OR THE INFORMATION, OR ANY DELAY OR INABILITY TO USE OR ACCESS THIS WEBSITE OR ANY OF THE INFORMATION, INCLUDING, WITHOUT LIMITATION, DIRECT, INDIRECT, INCIDENTAL, SPECIAL, CONSEQUENTIAL OR PUNITIVE DAMAGES, WHETHER UNDER A CONTRACT, TORT, STRICT LIABILITY OR ANY OTHER THEORY OF LIABILITY, EVEN IF THE COMPANY IS AWARE OR HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. THE COMPANY ALSO ASSUMES NO RESPONSIBILITY, AND SHALL NOT BE LIABLE FOR, ANY DAMAGES OF ANY KIND WHATSOEVER IN CONNECTION WITH OR RELATED TO ANY DAMAGE TO, OR VIRUSES THAT MAY INFECT, YOUR COMPUTER EQUIPMENT OR OTHER PROPERTY ON ACCOUNT OF YOUR ACCESS TO, USE OF, OR BROWSING OF THIS WEBSITE OR YOUR USE OR DOWNLOADING OF ANY INFORMATION FROM THIS WEBSITE.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

- TO THE EXTENT THIS WEBSITE IS BEING USED TO ACCESS ANY SOFTWARE OR APPLICATIONS PURSUANT TO RIGHTS GRANTED TO YOU UNDER OTHER AGREEMENTS ENTERED INTO BY YOU AND THE COMPANY, THEN YOUR USE OF THIS WEBSITE IS ALSO SUBJECT TO ALL DISCLAIMERS AND LIMITATIONS OF LIABILITY SET FORTH IN SUCH OTHER AGREEMENTS BETWEEN YOU AND THE COMPANY WHICH PROVIDE YOU ACCESS TO SUCH SOFTWARE OR APPLICATIONS. TO THE EXTENT OF ANY CONFLICT, AMBIGUITY OR INCONSISTENCY BETWEEN THESE TERMS OF USE AND THE TERMS OF ANY SUCH SPECIFIC OTHER AGREEMENTS, THE TERMS OF SUCH SPECIFIC OTHER AGREEMENTS SHALL CONTROL TO THE EXTENT OF SUCH CONFLICT, AMBIGUITY OR INCONSISTENCY.
3. No Unlawful Use. You may not use this Website for any purpose that is unlawful or prohibited by these terms and conditions.
 4. Links to Other Websites. This Website may, from time to time, contain hyperlinks to Websites operated or with content provided by third parties. Such other Websites are provided for your use and reference only and your use thereof is at your own risk. The Company does not control such Websites, and is not responsible for their contents. The inclusion of hyperlinks to such Websites and your use thereof does not imply any endorsement of the material on such Websites or any association with their operators on the part of the Company.
 5. Communications with The Company. In the event you transmit any communication or material to our Website by electronic mail or otherwise, including, without limitation, any data, questions, comments, suggestions or similar items, such transmission or communication is, and will be treated as, non-confidential and non-proprietary. Any such communications or materials you transmit or post on our Website becomes the property of the Company and may be used by the Company for any purpose, including, without limitation, reproduction, disclosure or transmission to other parties, publication, broadcast and posting. You are prohibited from posting or transmitting any unlawful, threatening, libelous, defamatory, obscene, scandalous, inflammatory, pornographic, or profane material or any material that could constitute or encourage conduct that would be considered a criminal offense, give rise to civil liability, or otherwise violate any law. The Company will cooperate with any law enforcement authorities or court order requesting or directing the Company to disclose the identity of anyone posting any such information or materials.
 6. Trademarks. The trademarks, logos, graphics and service marks (collectively the “Intellectual Property Rights”) displayed on our Website are registered and unregistered Intellectual Property Rights of the Company or its affiliates and third party vendors. Nothing contained on our Website grants or should be construed as granting, by implication, estoppel, or otherwise, any license or right to use any of the Intellectual Property Rights on our Website, without the express prior written consent of the Company or such affiliates and third party vendors that may own such Intellectual Property Rights, in each instance, which consents may be granted or withheld in the sole and absolute discretion of the Company. Except as specifically provided in these terms and conditions, your use of any of the Intellectual Property Rights, or any other information or content displayed on our Website is strictly prohibited. The Company will enforce its intellectual property rights to the fullest extent of the law, including the seeking of criminal prosecution.
 7. Terms and Conditions Subject to Change. The Company reserves the right to at any time revise these terms and conditions. Your continued use of this Website after any such revisions will constitute your acceptance of such revised terms and conditions. We encourage you to periodically visit this page to review the then current terms and conditions to which you are bound.
 8. General. These terms and conditions will be construed in accordance with the laws of the State of Florida, without regard to its conflict of law provisions. The provisions of these terms and conditions are severable. If any such provision is determined to be unenforceable, such unenforceability will not affect any other provision of these terms and conditions, and the terms and conditions will be construed as if such unenforceable provision had not been contained herein.

- No waiver or consent by the Company of any right under or provision of these terms and conditions will be deemed to be either a waiver of any other such right or provision or a waiver of that same right or provision at any other time.
9. Global Privacy Policy. OmniComm Systems respects and follows the best practices and privacy principles of Notice, Choice, Accountability for Onward Transfer, Security, Data Integrity and Purpose Limitation, Access, and Recourse, Enforcement and Liability including the transfer of personal information in electronic, paper or verbal formats from the member states of the European Union and Switzerland to the United States as set forth in the EU-U.S. Privacy Shield and the Swiss-U.S. Privacy Shield Framework.

<http://www.omnicomm.com/privacy>

To confirm participating organizations you may refer to <https://www.privacyshield.gov/list>

All content copyright © 2000-2018 OmniComm Systems, Inc. All rights reserved. For any information concerning the proprietary nature of the information contained on this website, please contact OmniComm's General Counsel at 954-473-1254 or e-mail legal@omnicomm.com.

2.0 Intended Audience

The prime users for AutoEncoder are Medical Coders. Other users include Data Managers and Monitors who require the ability to access, evaluate, and apply coding from industry dictionaries.


3.0 Manual Conventions

After each instruction, feedback is provided as necessary (i.e., if a window displays or a specific action occurs).

Bolded text highlights certain screenshots or indicates a user action, such as clicking selecting an option. In screenshots, user actions are outlined in red boxes or circles.

Example:

1. Click the **User Account** icon on the top right of the AutoEncoder status bar.



2. Select **Change Password** from the drop-down list.
3. Enter the **Old Password**, **New Password**. (Refer to [Password Rules](#) when entering the new password.)
4. Enter the new password in the **Confirm Password** field.
5. Click on the **Change Password** button to save the new password. A confirmation message displays at the top of the Login page and an email is sent to you, indicating that the password was successfully changed.

4.0 Getting Started


After logging into AutoEncoder refer to the [Autoencoder Workflow](#) to determine your tasks according to your user/role.

To get acquainted with the AutoEncoder interface, see [AutoEncoder Interface](#).

LOG IN/LOG OFF

4.1 Log In


The initial installation includes a single System User/System Admin (Administrator) named “AecAdmin” (not case-sensitive). Then you may then create other user accounts, which may be System Users/System Admins or Regular Users. (See [Edit Add a User](#) for more information.)

 **Note:** If you get locked out for exceeding the allowed number of incorrect login attempts, you must contact the trial administrator to have your password reset immediately. Otherwise, you cannot log in until an hour passes.

To log in:


1. Start a web browser application, such as IE 11, Firefox, or Chrome.
2. Enter the appropriate URL for AutoEncoder. The URL is set up when AutoEncoder is first installed.

The Log in page displays.



The screenshot shows the 'AutoEncoder Login' page. It has a title 'AutoEncoder Login' in blue. Below it is a label 'User Name or Email' followed by a text input field with the placeholder 'Username'. Below that is a label 'Password' followed by a text input field with the placeholder 'Password'. At the bottom, there is a blue 'Login' button and a link 'Forgot Password?' in blue text.

3. Enter your User Name/Email address and Password and click the **Login** button.
4. If you were assigned multiple roles, click on a role in the Select Role page.



The screenshot shows the 'AutoEncoder: Select Role' page. It has a title 'AutoEncoder: Select Role' in blue. Below the title, there is a list of roles, each preceded by a blue circular icon with a right-pointing arrow. The roles are: Read, Coder, Reviewer, Approver, Super User, and ReadOnly Monitor.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

The Dashboard displays.

AutoEncoder Dashboard

Select a module:

▲ Trial	Module	Dicti...	Code	Review	Approve	Organizat...	Trial Group	Actions	Job Status
Demo Study 1	WHODRG	WhoDrug...	2	0	0	OmniComm		Reports	100% Completed
Demo Study 1	MEDRA	MedDRA ...	20	1	0	OmniComm		Reports	
Demo Study 2	WHODRG	WhoDrug...	0	0	0	OmniComm		Reports	
Demo Study 2	MEDRA	MedDRA ...	9	0	0	OmniComm		Reports	
Demo Study Edi...	MEDRA	MedDRA ...	0	-	-	OmniComm	FLA	Reports	
P500-0717	WHODRG	WhoDrug...	0	-	-	FLA	TG A	Reports	
SDTM_DEMO	MH	MedDRA ...	9	1	1	OmniComm	FLA	Reports	
SDTM_DEMO	CPCODE	MedDRA ...	1	0	0	OmniComm	FLA	Reports	
SDTM_DEMO	MedDRA	MedDRA ...	9	1	0	OmniComm	FLA	Reports	
SDTM_DEMO	WHODRG	WhoDrug...	6	0	4	OmniComm	FLA	Reports	100% Completed
TrialA	AECODE	MedDRA ...	98	0	0	Organization A	TG A	Reports	

14 items

Note: If the user role has no access to the data on the page, no rows/data display. However when navigating from the Coding and Transmission Queue pages, the Dashboard displays the relevant data.

5. Refer to [Use the Dashboard for Coding](#) to get started with AutoEncoder.

4.2 Change Password


For security purposes, it may be necessary to change your password. There are two ways to change your password:

- Change your password manually through the User Account Manager
- Wait until the system prompts you for a new password when it expires. Once your password expiration is near, a warning displays after login. The system gives you the option to change it immediately, or sometime before the password expiration (the default setting is 5 days, which may be configured to another setting).

Note: If your password becomes locked, reset your password using the **Forgot Password** link on the login page. (See [Retrieve a Forgotten Password](#) for instructions.)


To change your password:

1. Click the **User Account** icon on the top right of the AutoEncoder status bar.




2. Select **Change Password** from the drop-down list.

The Change Password dialog displays.



3. Enter the **Old Password** and **New Password**. (Refer to [Password Rules](#) when entering the new password.)
4. Enter the new password in the **Confirm Password** field.
5. Click on the **Change Password** button to save the new password. A confirmation message displays at the top of the Login page and an email is sent to you, indicating that the password was successfully changed.
6. Click the **Cancel** button to return to the Login page.

 **Note:** If you do not use the link your password remains the same. If you did not request your password be reset, please contact your administrator.

4.2.1 Password Rules

The following password rules apply (these settings are configurable):


- Must have at least 1 uppercase letter
- Must have at least 1 lowercase letter
- Must have at least 1 numeric character
- Password length is must be 8 characters
- Must have at least 1 special character
- Password expiration is 90 days
- Password Warning message displays on the Dashboard on login starting 5 days prior to the expiration:
"Your password will expire in {5} days. Please change it as soon as possible."

4.2.2 Retrieve a Forgotten Password


When retrieving a forgotten password, a temporary password is sent to your email address.

To retrieve a forgotten password:

1. Click the **Forgot Password** link on the Login page.

The image shows the 'AutoEncoder Login' dialog box. It has a title bar 'AutoEncoder Login'. Below the title bar, there are two input fields: 'User Name or Email' with the text 'lvioleta' and 'Password' with the text 'Password'. At the bottom, there are two buttons: 'Login' and 'Forgot Password?'. The 'Forgot Password?' button is circled in red.

The Recover Password dialog displays.

The image shows the 'Recover Password' dialog box. It has a title bar 'Recover Password'. Below the title bar, there is a label 'Email address:' followed by an input field containing the text 'lvioleta@omnicomm.com'. At the bottom, there are two buttons: 'Recover Password' and 'Cancel'. The 'Recover Password' button is highlighted with a red border.

2. Enter your email address and click the **Recover Password** button. The Login page re-displays and you receive two notification emails: A notification with a reset password link and a notification with a temporary password.
3. Enter your **User Name** and the temporary **Password** in the Login page. (The password can only be used once and you must change it as soon as possible.)
4. Enter the required information in the Change Password dialog. (Refer to [Change Password](#) for guidance.)

4.2.3 Password Notifications

Temporary password notifications are sent when a:

- New user is created
- System User/System Admin clicks the Reset Password link on the Users page to reset a user's password
- User clicks the **Forgot Password** link on the Login page

The system sends an email notification with a newly generated temporary password to the user. The password is valid for the configured password expiration (default is 90 days). This password is a one-time use password and must be changed after the user logs in.

Users are not forced to change the password on login. If the user does not change the password after login, they will not be able to log in again.


The user can click on the **Forgot Password** link on the login page or the System User/System Admin can click on the **Reset Password** link for that user to reset the user's password.

4.2.4 Log Off

You can log off from the User Account drop-down menu.

To log off:

1. Click the **User Account** icon on the top right of the AutoEncoder status bar.



2. Select **Log Off** from the drop-down list.

USING AUTOENCODER

5.0 About AutoEncoder

Most data entered into EDC (Electronic Data Capture) is entered as free text but must be coded against standard medical dictionaries. Commonly coded clinical data includes AE (Adverse Event) verbatim terms, medical history, physical conditions and procedures, and medications.

To start using AutoEncoder, refer to [AutoEncoder Workflow](#).

6.0 AutoEncoder Workflow


This section provides AutoEncoder workflow diagram and steps.

AutoEncoder workflow steps:


1. The System Administrator does the following:
 - a. Creates needed objects such as organizations, trial groups, and trials (see [Manage Organizations](#), [Manage Trial Groups](#), or [Manage Trials](#))
 - b. Loads all needed Dictionary Versions (see [Add New Dictionary Version](#))
 - c. Manages System-Wide Settings (see [View/Edit System Settings](#))
 - d. Bulk imports any required Synonym Lists (see [Import/Export Synonym List](#))
 - e. Creates any needed users (see [Edit/Add a User](#))
 - f. Assigns user(s) to appropriate objects with a corresponding role (see [Assign a User to an Object](#))
2. The System Administrator may assign a user the role "Super User" on an organization (see [Assign User to an Object](#)). The Super User, when assigned to an organization, does the following in the context of their parent organization:
 - a. Creates new users (see [Edit/Add a User](#))
 - b. Assigns user(s) to any object belonging to the Super User's organization (see [Assign User to an Object](#))
 - c. Edits users (see [Edit/Add a User](#))
 - d. Assign user(s) to any object owned by the organization as any applicable role (see [Assign a User to an Object](#))
 - e. Associates dictionary versions with modules (see [Manage Mapping](#))
 - f. Matches/maps TrialMaster module and AutoEncoder dictionary terms (see [Manage Mapping](#))
 - g. Manages Synonym Lists (see [Manage Synonym Lists](#))
 - h. Manages settings on individual objects (see [Manage Organizations](#), [Manage Trial Groups](#), or [Manage Trials](#))
3. Once a module and its associated mapping have been established, the process of autocoding and manual validation begins. Regular users with appropriate permissions can:
 - View, confirm, review and/or approve terms after they load into AutoEncoder (see [Confirm/Review/Approve Terms](#))
 - Manually code and uncode existing terms (see [Code/Unicode Matched Terms](#))
 - Manually uncode pre-filled terms already passed back to the associated application (see [Code/Unicode Matched Terms](#)).

4. At any point, Regular Users with the appropriate permission can view reports and view the audit history for individual terms (or various management/objects). (See [View Module Reports](#), [View System Reports](#), and [View a Term's Audit History](#).)

The System User/System Admin, Super User and Regular User workflow diagrams illustrate the typical workflow process for each user type.



AE System User/System Admin Workflow



7.0 AutoEncoder Interface

The AutoEncoder window consists of the following:

- [Dashboard \(Home Page\)](#)
- [Menus](#)
- [Status Bar](#)

Main AutoEncoder Window

7.1 Home Page (Dashboard)

Once you log in the Dashboard displays. See [Using the Dashboard for Coding](#) for instructions.

AutoEncoder Dashboard


Select a module: 🔄 ✎ 📄 🔍

▲ Trial	Module	Dicti...	Code	Review	Approve	Organizat...	Trial Group	Actions	Job Status
Demo Study 1	WHODRG	WhoDrug...	2	0	0	OmniComm		Reports	100% Completed
Demo Study 1	MEDRA	MedDRA ...	20	1	0	OmniComm		Reports	
Demo Study 2	WHODRG	WhoDrug...	0	0	0	OmniComm		Reports	
Demo Study 2	MEDRA	MedDRA ...	9	0	0	OmniComm		Reports	
Demo Study Edi...	MEDRA	MedDRA ...	0	-	-	OmniComm	FLA	Reports	
P500-0717	WHODRG	WhoDrug...	0	-	-	FLA	TG A	Reports	
SDTM_DEMO	MH	MedDRA ...	9	1	1	OmniComm	FLA	Reports	
SDTM_DEMO	CPCODE	MedDRA ...	1	0	0	OmniComm	FLA	Reports	
SDTM_DEMO	MedDRA	MedDRA ...	9	1	0	OmniComm	FLA	Reports	
SDTM_DEMO	WHODRG	WhoDrug...	6	0	4	OmniComm	FLA	Reports	100% Completed
TrialA	AECODE	MedDRA ...	98	0	0	Organization A	TG A	Reports	

14 items

7.2 Status Bar

The Status Bar displays the following information (from left to right):




- Main Menu icon – provides access to the main sub menus
- AutoEncoder logo – allows you to return to the home page (Dashboard)
- Active/Timeout – Either displays the time remaining in the session (if the session is inactive) or an “Active” icon when the session is in use
- User – identifies the user name currently logged in
- Role – identifies the role of the user
- User Account icon – allows you to edit your account information, change your password, or log off


7.3 Menus

Note: The menus displayed vary depending upon permissions assigned to the active user and role.


The following menus display in Autoencoder:

- Main Menu – provides access to the main sub menus



Click the  **Main Menu** icon (top left) to access specific pages such as Dashboard, Trial Groups, Organizations, and Users.

- User Account




Click the  **User Account** icon to do the following:


- Edit your account information
- Change your password
- Log off

Note: See [Update Your User Account Information](#) for more information about your user account.

- Information Menu (bottom of main menu) – provides copyright notice and build number and the following links:



- **Help** – displays the AutoEncoder Help window
- **About** – displays the **About AutoEncoder page**, which provides the system message, environment, copyright, and database and software versions



- **Contact** – displays the OmniComm Customer Care page, which provides contact information for different customer regions. It also includes the customer care email

Contact OmniComm Customer Care

OmniComm Customer Care

Contact information:

North America Toll Free: (866) 996-6332

International Toll Free: +800-OMNICOMM

International Toll Direct: +49-228-227-44-11

India Toll Free: 000800-1006778

Mexico Toll Free: 001800-05140455

Panama Toll Free: 00800-2265719

Ukraine Toll Free: 0800-5041810

(Specific toll free numbers can be arranged if required)

Email: customercare@omnicomm.com

- **Terms** – displays the standard user agreement


Dashboard

7.4 Use the Dashboard for Coding

The Dashboard page displays the dictionary, number of terms that need to be coded/reviewed/approved, organization, trial group, and job status for each module.

The following links are provided to view more detailed information:

- [Module/Code/Review/Approve](#) – links to the Coding page of the applicable module where you can [confirm/review/approve terms](#)
- [Actions](#) – links to reports of the applicable module
- [Job Status](#) – links to the Autoencoding Jobs page where you can view more details about each autoencoding job

 **Note:** If a “-” appears in the Review or Approve column rather than a numerical value, this means that AutoReviewApprove (Automatically Review and Approve Terms once Confirmed) setting is on (set to True), and therefore, all items are accessible from the coding stage (Code column). The only exception is when AutoReviewApprove setting is turned on mid-stream. In this case, there could be remnant terms in the Review and/or Approve Stage and once these terms complete the cycle the “-” appears. (See [View/Edit System Settings](#) for more information.)

The Dashboard allows you to do the following, depending on your role permissions:

- [View coding stages](#)
- [View a term's audit history](#)
- [Confirm/review/approve terms](#)
- [Re-autoencode terms](#)
- [Code/decode matched terms](#)
- [View AutoEncoding job details](#)
- [View/decode terms in the Transmission Queue](#)
- [View Module Reports](#)

7.4.1 View Coding Stages

You can view coding stages (coding, review, approve) for each group of terms in each module.

To view coding stages:

1. Access the Dashboard after first logging in or by clicking the AutoEncoder logo.



The Dashboard page displays.

AutoEncoder Dashboard

Select a module: ↺ ↻ 🔍 📄

▲ Trial	Module	Dicti...	Code	Review	Approve	Organizat...	Trial Group	Actions	Job Status
Demo Study 1	WHODRG	WhoDrug...	2	0	0	OmniComm		Reports	100% Completed
Demo Study 1	MEDRA	MedDRA ...	20	1	0	OmniComm		Reports	
Demo Study 2	WHODRG	WhoDrug...	0	0	0	OmniComm		Reports	
Demo Study 2	MEDRA	MedDRA ...	9	0	0	OmniComm		Reports	
Demo Study Edi...	MEDRA	MedDRA ...	0	-	-	OmniComm	FLA	Reports	
P500-0717	WHODRG	WhoDrug...	0	-	-	FLA	TG A	Reports	
SDTM_DEMO	MH	MedDRA ...	9	1	1	OmniComm	FLA	Reports	
SDTM_DEMO	CPCODE	MedDRA ...	1	0	0	OmniComm	FLA	Reports	
SDTM_DEMO	MedDRA	MedDRA ...	9	1	0	OmniComm	FLA	Reports	
SDTM_DEMO	WHODRG	WhoDrug...	6	0	4	OmniComm	FLA	Reports	100% Completed
TrialA	AECDCE	MedDRA ...	98	0	0	Organization A	TG A	Reports	

14 items

2. Do one of the following:
 - To access all stages, click the applicable link under the **Module** column.
 - To access a specific stage, click the applicable link under the **Code**, **Approve**, or **Review** column. (The Approve and Review links display a "-" if the AutoReviewApprove setting is set to True.)

The Coding page displays the verbatim term and any corresponding reference items with each patient, site, visit, and form. (These first four columns are specific to TrialMaster.) Links are provided to view the associated TrialMaster form, match (to manual/autocode/history/synonym/ dictionary), and the date the term was last modified (audit history).

Coding: SDTM_DEMO - MH Transmission Queue

☒ Approve Stage 1 item
Approve



Patient	Site	Visit	Form	MedDRA LLT Name	▲ Match	LLT Name	LLT Cc
<input type="checkbox"/> 100-101-...	100	SCREEN	MH	High LDL	Manual - Dictionary Match	Muzzy head	10028


☒ Review Stage 1 item
Review Approve


Patient	Site	Visit	Form	MedDRA LLT Name	▲ Match	LLT Name	LLT Cc
<input type="checkbox"/> 123-456-...	123	SCREEN	MH	Pyoderma Gangrenosum	Manual - Dictionary Match		

☒ Coding Stage 9 items Re-autocode
Confirm Review Approve

Patient	Site	Visit	Form	MedDRA LLT Name	▲ Match	LLT Name	LLT Cc
<input type="checkbox"/> 100-101-...	100	SCREEN	MH	Appendectomy	Manual - Dictionary Match	Cauliflower ear	10
<input type="checkbox"/> 123-678-...	123	SCREEN	MH	Pyoderma Gangrenosum	Manual - Dictionary Match	Pyogenic sterile arthritis pyo...	10
<input type="checkbox"/> 123-456-...	123	SCREEN	MH	Space Test	No Match		
<input type="checkbox"/> 123-456-...	123	SCREEN	MH	Space Test	No Match		

Each stage is automatically expanded so all data is viewable. The header of each stage is highlighted with a color and indicates the number of terms included in the stage. You can expand/collapse a stage using the  Expand/ Collapse icon.

 **Note:** If a “-” displays in the Review or Approve column rather than a numerical value means that AutoReviewApprove (*Automatically Review and Approve Terms once Confirmed*) setting is On, and therefore, all items are accessible from the Coding Stage (Code column). The only exception is when AutoReviewApprove setting is turned on mid-stream. In this case, there could be remnant terms in the Review and/or Approve Stage and once these terms complete the cycle the “-” displays. (See [View/Edit System Settings](#) for more information.)

 **Note:** You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions](#).)

7.4.2 View a Term's Audit History

You can view a term's audit history from the following locations:

- Coding, Review, or Approve stage
- Transmission Queue

To view a term's audit history:

1. Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
2. Do one of the following:
 - Locate the applicable term in the Coding, Review, or Approve stage.
 - Click the **Transmission Queue** button and locate the term.
3. Click the link of the term under the **Last Modified** column (far right).

Coding: Trial B - AE Transmission Queue

Coding Stage 10 Items Re-autocode Confirm Review Approve ↺ ↻ ✕ ☰					
LLT Name	LLT Code	PT Name	PT Code	Last Modified	
Headache dull	10019215	Headache	10019211	27-NOV-2017 07:16:03 AM	
Headache	10019211	Headache	10019211	14-NOV-2017 01:15:38 AM	
Chronic gastric ulcer without...	10055768	Gastric ulcer, obstructive	10017840	27-NOV-2017 02:02:50 AM	
Headache	10019211	Headache	10019211	14-NOV-2017 01:19:58 AM	

The Term Audit Report displays the audit history of the term, including the user who performed the action, the trial, module name, status, match, and dictionary term. (If the System Setting, "IncludeTimeZoneInfo" is set to True, the timezone information is appended to all of the timestamps. (See [View/Edit System Settings](#) for more information.)

Term Audit Report Back

Term Audit Report ↺ ↻ ✕ ☰								
Audit Date	By User	Trial Name	Trial Description	Module Name	Term	Status	Match	
27-NOV-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Approved	Manua	▲
14-NOV-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Reviewed	Manua	
18-OCT-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Confirmed	Manua	
17-OCT-2017 0...	Marcela F(49)	Trial B	Trial B2	AE	Dull headache	Assigned	Manua	
13-OCT-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Confirmed	Manua	
13-OCT-2017 1...	Marcela F(49)	Trial B	Trial B2	AE	Dull headache	Assigned	Manua	
13-OCT-2017 1...	Marcela F(49)	Trial B	Trial B2	AE	Dull headache	Confirmed	Manua	
13-OCT-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Assigned	Manua	
22-AUG-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Assigned	Autoco	
17-MAY-2017 0...	AutoEncoder A...	Trial B	Trial B2	AE	Dull headache	Assigned	Autoco	

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.


Note: You can sort and filter the report or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))

7.4.3 Confirm/Review/Approve Terms

The following actions allow you to move a term to the next Coding Stage:

- Confirm – moves term to Review Stage
- Review – moves item to Approve Stage
- Approve – moves item to Transmit Stage

The Coding Stages diagram illustrates how to move to each stage. As shown in this diagram, you can also re-autocode terms from the Coding Stage before moving them to the Review stage.



AutoEncoder is configured so you can skip the Review or Approve Stage via the AutoReviewApprove setting. However, you can change this setting to go through each stage (see [Edit Object Settings](#) and [View/Edit System Settings](#)).

To confirm/review/approve terms:

1. Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
2. Click the check boxes of the applicable terms under the stage header. (To select all of the items, click the **Select All** check box just below the table header.)

Coding: Trial A - AE Transmission Queue

☒ **Approve Stage** 2 Items Approve

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Code
<input type="checkbox"/> 9933	Miami C...	This Visit	AE #1	Gangrene	<input checked="" type="checkbox"/> Manual - Dictionary Match	Jejunal gangrene	10023172
<input type="checkbox"/> 5623	Texas Te...	Other Visit	AE #1	Pain in head	<input checked="" type="checkbox"/> Manual - History Match	Headache dull	10019215

☒ **Review Stage** 2 Items Review Approve

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Code
<input type="checkbox"/> 9933	Miami C...	This Visit	AE #1	Gangrene	<input checked="" type="checkbox"/> Autocode - Dictionary Match	Gangrene	10017771
<input type="checkbox"/> 2121	Miami C...	Other Visit	AE #1	Dizziness	<input checked="" type="checkbox"/> Manual - Dictionary Match	Headache dull	10019215


☒ **Coding Stage** 8 Items Re-autocode **Confirm (2)** Review (2) Approve (2)

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Code
<input checked="" type="checkbox"/> 2344	Miami C...	That Visit	AE #1	Broken arm	<input checked="" type="checkbox"/> Manual - Dictionary Match	Aspirin asthma	10003172
<input checked="" type="checkbox"/> 3342	Texas Te...	Other Visit	AE #1	Dull headache	<input checked="" type="checkbox"/> Manual - Dictionary Match	Headache dull	10019215

The terms and the applicable button in the header (Confirm/Review/Approve) become highlighted. The header indicates the number of terms you selected. For example, the Confirm button displays when selecting items under the Coding Stage header, as shown in the figure above.

3. Select the **Confirm/Review/Approve** button. (From the Coding Stage, you can click the **Approve** button to skip the Review Stage.)

A confirmation message displays at the top of the page and the term is moved to the next stage. For example, if you confirmed the terms, the terms move to the Review Stage. However, if you selected to approve the terms, they are moved to the Transmission Queue.



Coding: Trial A - AE [Transmission Queue]

Approve Stage 3 Items [Approve] [Refresh] [Print] [Export]

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Co
6588	Texas Te...	Other Visit	AE #2	Headake	Manual - Dictionary Match	Headache occurring	10
9933	Miami C...	This Visit	AE #1	Gangrene	Manual - Dictionary Match	Jejunal gangrene	10

Review Stage 3 Items [Review] [Refresh] [Print] [Export]

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Co
2344	Miami C...	That Visit	AE #1	Broken arm	Manual - Dictionary Match	Aspirin asthma	10
3342	Texas Te...	Other Visit	AE #1	Dull headache	Manual - Dictionary Match	Headache dull	10

Coding Stage 6 Items [Re-autocode] [Confirm] [Review] [Approve] [Refresh] [Print] [Export]

Patient	Site	Visit	Form	MedDRA LLTName	Match	LLT Name	LLT Co
7272	Texas Te...	That Visit	AE #3	Persistent headache	Manual - Dictionary Match	Aged in-law	10
9378	Texas Te...	Other Visit	AE #3	ear canal pimple	Manual - Dictionary Match	Ear swelling	10

- If you need manually uncode/code a term, see [Code/Unicode Matched Terms](#).


7.4.4 Re-autocode Terms

Re-autocoding is only necessary if you added new synonyms or history terms and then you want to re-code no-match items against them. Normally, autocoding occurs automatically against a term when it is pulled from another application such as TrialMaster.

Note: The re-autocode feature matches the entire module.

To re-autocode terms:

- Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
- Click the **Re-autocode** button under the Coding stage.



Coding Stage 28 Items [Re-autocode] [Confirm] [Review] [Approve] [Refresh] [Print] [Export]

Patient	Site	Visit	Form	Trade Name	Route for Coding	Indication for Coding	Co
123-678...	123	CON ME...	WHODRG	Depakene	SUBCUTANEOUS	PROPHYLAXIS/SUPPLEMENT	
123-456...	123	CON ME...	WHODRG	Vancomycin	TOPICAL	Pyoderma Gangrenosum	
123-678...	123	CON ME...	WHODRG	Pharnolone	INTRAMUSCULAR	PROPHYLAXIS/SUPPLEMENT	
100-111...	100	CON ME...	WHODRG	BUSONAL	RESPIRATORY (INHALATION)	PROPHYLAXIS/SUPPLEMENT	

A green success message displays at the top of the page, indicating the job number and that the autocoding job was started. A red failure message displays at the top of the page if there is a job already in the queue.

- To check the progress of the job, return to the Dashboard and locate the applicable link under the Job Status column on the left-hand side of the table. The link displays the progress percentage between 0 and 100. (Click the link to view more details.)

AutoEncoder Dashboard

Select a module:

▲ Trial	Module	Dictionary	Code	Review	Approve	Organization	Trial Group	Actions	Job Status
Trial A	AE	MedDRA 16.1	0	4	3	Customer A	Trials A and B	Reports	100% Completed
Trial A	CONMED	WhoDrug B2 Sep 2016	0	1	1	Customer A	Trials A and B	Reports	100% Completed
Trial B	AE	MedDRA 16.1	8	1	1	Customer B	Trials A and B	Reports	100% Completed
Trial B	CONMED	WhoDrug B2 Sep 2013	8	0	1	Customer B	Trials A and B	Reports	100% Completed

7.4.5 Code/Unicode Matched Terms

You can code/unicode the following types of matched terms:

- No Match
- Autocode – History Match
- Autocode – Synonym Match
- Autocode – Dictionary Match
- Manual – History Match
- Manual – Synonym Match
- Manual – Dictionary Match

To code/unicode matched terms:

- Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
- Click the applicable link under the **Match** column.

Coding: Trial B - MedDRA

Approve Stage: 3 Items

Patient	Site	Visit	Form	MedDRA LLT Name	Match	LLT Name	LLT Code	SOC Name	PT Name
003-006-CRB-005	003	Adverse	MEDDRA	Headache	Autocode - Dictionary Match	Headache	10019211	Nervous system disorders	Headache
003-008-JER-007	003	Adverse	MEDDRA	Fatigue	Manual - Dictionary Match	Fatigue management	10016258	Surgical and medical proced...	Fatigue man...

The header of the Manual Coding page displays the verbatim term plus any reference values and application column values. The middle section displays the term that is currently selected and being coded, if there was a dictionary match, or it was autocoded. The bottom section displays the Dictionary Search table where you can search for another verbatim term you want to code/unicode.

Manual Coding: Trial B - MedDRA

Verbatim Term	Patient	Site	Visit	Form
Headache	003-006-CRB-005	003	Adverse	MEDDRA

Currently Selected Coded by AecAdmin at 11-DEC-2017 07:39:01 AM (UTC-05:00) Eastern Time (US & Canada)

Type	LLT Name	LLT Code	PT Name	PT Code	HLT Name	HLT Code	HLGT Name	HLGT Code	SOC Name	SOC Code
Dictionary	Headache	10019211	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200

Code **Uncode** Cancel

Dictionary Search 67 item matches: Verbatim Term like 'Headache' **Clear**

Type	Headache	LLT Code	PT Name	PT Code	HLT Name	HLT Code	HLGT Name	HLGT Code	SOC Name	SOC Code
Dictionary	Headache	10019211	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Throbbing headache	10058140	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Lumbar puncture headache	10025001	Post lumbar pu...	10060854	Neurological an...	10029294	Procedural relat...	10069888	Injury, poisonin...	10029200
Dictionary	Cervicogenic headache	10064888	Cervicogenic h...	10064888	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Headache drug withdrawal	10019214	Drug withdrawa...	10013753	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200

Note: If you clicked a **No Match** link, the middle section is empty and you can search for a term in the bottom section.

- If you want to search for another verbatim term, type the applicable search parameters in the corresponding fields and tab out or press the ENTER key. For example, to search for a migraine headache that is a nervous system disorder, type "Migraine Headache" in the first field (LLT Name) and "nervous system disorders" in the SOC Name field.

Manual Coding: Trial B - MedDRA

Verbatim Term	Patient	Site	Visit	Form
Headache	003-006-CRB-005	003	Adverse	MEDDRA

Currently Selected Coded by AecAdmin at 11-DEC-2017 07:39:00 AM (UTC-05:00) Eastern Time (US & Canada)

Type	LLT Name	LLT Code	PT Name	PT Code	HLT Name	HLT Code	HLGT Name	HLGT Code	SOC Name	SOC Code
Dictionary	Headache	10019211	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200

Code **Uncode** Cancel

Dictionary Search 5 items match: LLT Name like 'Migraine Headache', SOC Name like 'nervous system' **Clear**

Type	Migraine Headache	LLT Code	PT Name	PT Code	HLT Name	HLT Code	HLGT Name	HLGT Code	nervous syste...	SOC Code
Dictionary	Migraine headache	10027602	Migraine	10027599	Migraine heada...	10027603	Headaches	10019231	Nervous syste...	10029200
Dictionary	Headache (excl migraine) ag...	10056791	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Tension headache (excl migr...	10043270	Tension headac...	10043269	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Headache (except migraine)...	10019212	Headache	10019211	Headaches NEC	10019233	Headaches	10019231	Nervous syste...	10029200
Dictionary	Migraine type headaches	10027606	Migraine	10027599	Migraine heada...	10027603	Headaches	10019231	Nervous syste...	10029200

Note: By default, the dictionary search does an auto search using a like function as indicated in the header of the bottom section. You can clear the search by clicking on the **Clear** button or revert back to the original selection by clicking the **Revert** button (top right of page). Mapped reference columns are pre-filled in the search area.

- With the term selected (highlighted), click the **Uncode/Code** button. If you want to uncode a searched term in the Dictionary Search table, select (highlight) the term and then click the **Uncode** button. The Coding page re-displays and the term you uncoded/coded displays under the Coding Stage. In the example below, the LLT Name was sorted in ascending order by clicking the column header, and the SOC Name was dragged to the left so it is viewable without scrolling. (See [Sort and Filter Tables](#) for more information on sorting/filtering tables.)

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

Coding Stage 790 Items Re-autocode									
Patient	Site	Visit	Form	MedDRA LLT Name	Match	LLT Name	LLT Code	SOC Name	PT Name
003-006...	003	Adverse	MEDDRA	Headache	No Match				
003-006...	003	Adverse	MEDDRA	Headache	Manual - Dictionary Match	Migraine headache	10027602	Nervous system disorders	Migraine
003-009...	003	Adverse	MEDDRA	Headache	Autocode - Dictionary Match	Headache	10019211	Nervous system disorders	Headache
003-003...	003	Adverse	MEDDRA	Headache	Autocode - Dictionary Match	Headache	10019211	Nervous system disorders	Headache
003-006...	003	Adverse	MEDDRA	Headache-intermittent	No Match				
003-001...	003	Adverse	MEDDRA	Headache-intermittent	No Match				
004-002...	004	Adverse	MEDDRA	Hematuria	Autocode - Dictionary Match	Hematuria	10019450	Renal and urinary disorders	Haematuria

Note: If the term is coded using a history match, an expander shows/hides the history origin. When the history origin is expanded/shown, a message displays, showing where the history match originated. The history term originates from <object1> History, inherited from <object2>. For example, the history term originates from TrialA History, inherited from trial.

7.4.6 View Autoencoding Job Details

You can view the details of autoencoding jobs. The status of each job may include:

- Submitted
- Processing
- Failed
- Finished

To view job status details:

- Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
- Click the applicable link under the Job Status column. The AutoEncoding Jobs page displays the number jobs and details about each job, including progress (between 0 and 100%), status, trial and module name, and date submitted.

AutoEncoding Jobs for Trial A/AE Back								
AutoEncoding Jobs for Trial A/AE 20 Items								
Job Id	Status	Status Date	Progress	Trial Name	Module Name	Description	User	Date ...
50	Finished	05-DEC-2017 09:54:03 AM...	100	Trial A	AE	Trial A	Ivioleta	05-DEC-2...
49	Finished	01-DEC-2017 04:22:03 PM...	100	Trial A	AE	Trial A	Ivioleta	01-DEC-2...
48	Finished	01-DEC-2017 02:22:03 PM...	100	Trial A	AE	Trial A	Ivioleta	01-DEC-2...
47	Finished	01-DEC-2017 02:16:03 PM...	100	Trial A	AE	Trial A	Ivioleta	01-DEC-2...
46	Finished	01-DEC-2017 10:56:03 AM...	100	Trial A	AE	Trial A	AecAdmin	01-DEC-2...
45	Finished	01-DEC-2017 10:54:06 AM...	100	Trial A	AE	Trial A	AecAdmin	01-DEC-2...
44	Finished	30-NOV-2017 12:32:03 PM...	100	Trial A	AE	Trial A	Ivioleta	30-NOV-2...
40	Finished	16-NOV-2017 01:30:17 PM...	100	Trial A	AE	Trial A	AecAdmin	16-NOV-2...
36	Finished	14-NOV-2017 01:10:13 AM...	100	Trial A	AE	Trial A	AecAdmin	14-NOV-2...
35	Finished	14-NOV-2017 01:08:10 AM...	100	Trial A	AE	Trial A	AecAdmin	14-NOV-2...
31	Finished	10-NOV-2017 07:58:14 AM...	100	Trial A	AE	Trial A	AecAdmin	10-NOV-2...

Note: You can sort and filter the table or export it to Excel. Various table options are also available.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

(See [Use Table Functions](#).)

7.4.7 View/Unicode Terms in Transmission Queue

The transmissions queue shows the transmit status, which can be one of the following

- Ready To Send – The term has reached Approved Status and is ready to be sent to the calling application (TrialMaster). This occurs automatically by the service.
- Sent – AutoEncoder service has sent it to the calling application and is awaiting a response.
 - If the transmission fails, it automatically retries for approximately 1 week: 1 time per hour for 150 times.
 - If there are no issues with the connection, the status changes from Ready to Send directly to Completed.
- Send Failed – AutoEncoder tried to send, but had an error and will retry
- Completed – AutoEncoder got confirmation that the calling application (TrialMaster) has been updated. At this point it is available for uncoding if desired.

To view the transmission queue:

1. Access the Dashboard and click a code, module, review, or approve link. (See [View Coding Stages](#) for guidance.)
2. Click the **Transmission Queue** button in the Coding page.

Coding: Trial B - MedDRA

Transmission Queue

Approve Stage3 items

<input type="checkbox"/>	Patient	Site	Visit	Form	MedDRA LLT Name	Match
<input type="checkbox"/>	001-001-...	001	Adverse	MEDDRA	Insomnia	Autocode - Dictionary Match
<input type="checkbox"/>	001-001-...	001	Adverse	MEDDRA	Anxiety	Autocode - Dictionary Match

Review Stage3 items

<input type="checkbox"/>	Patient	Site	Visit	Form	MedDRA LLT Name	Match
<input type="checkbox"/>	003-001-...	003	Adverse	MEDDRA	Febrile neutropenia	Autocode - Dictionary Match
<input type="checkbox"/>	003-001-...	003	Adverse	MEDDRA	Alanine aminotransferase in...	Autocode - Dictionary Match

The Transmission Queue page displays the transmit status, patient name, site, visit, form, and dictionary match for each term. (The Form column contains links to the associated forms.)

Transmission Queue: Trial B - MedDRA [Back](#)

Transmission Queue 5 items [Uncode \(2\)](#)

<input type="checkbox"/>	Transmit Status	Patient	Site	Visit	Form	MedDRA LLT Name	Match	LLT Name
<input type="checkbox"/>	Ready To Send	003-001-...	003	Adverse	MEDDRA	Dehydration	Autocode - Dictionary Match	Dehydration
<input checked="" type="checkbox"/>	Ready To Send	003-001-...	003	Adverse	MEDDRA	Anorexia	Autocode - Dictionary Match	Anorexia
<input checked="" type="checkbox"/>	Ready To Send	003-001-...	003	Adverse	MEDDRA	Alanine aminotransferase in...	Autocode - Dictionary Match	Alanine aminotransfer
<input type="checkbox"/>	Ready To Send	003-001-...	003	Adverse	MEDDRA	Hypoalbuminemia	Autocode - Dictionary Match	Hypoalbuminemia
<input type="checkbox"/>	Ready To Send	003-001-...	003	Adverse	MEDDRA	Sinus tachycardia	Autocode - Dictionary Match	Sinus tachycardia

- If you need to uncode terms, select the check box of the applicable terms and click the **Uncode** button (top right). (When you select the terms, the Uncode button becomes highlighted and displays the number of terms you selected.)

The number of terms selected for uncoding display at the top of the page.

- Click the **Back** button at the top of the page to return to the Coding page.

7.4.8 View Module Reports

You view the following module reports from the Dashboard:

- Coding
- Coding Transmission History
- Dictionary Mapping History
- Dictionary Version History
- Mapping

To view module reports:

- Access the AutoEncoder Dashboard after first logging in or by clicking the AutoEncoder logo.



- Click the **Reports** link of the applicable module under the Actions column (far right).

AutoEncoder Dashboard


Select a module:

Trial	Module	Dictionary	Code	Review	Approve	Organization	Trial Group	Actions	Job Status
Demo Study	MEDRA	MedDRA MedDRA ...	0	0	0	Healthy Solutions	Allergy	Reports	
Trial B	MedDRA	MedDRA MedDRA ...	795	3	3	OmniComm	Trial Group AB	Reports	100% Completed
Trial B	ST2	WhoDrug B3 Who ...	8	0	0	OmniComm	Trial Group AB	Reports	100% Completed

The Module Reports page lists module report links.



3. Click a link to review a report.

 **Note:** You can sort and filter the report or export it to Excel. Various table options are also available. (See [Use Table Functions](#).)

Basic AutoEncoder Functions


7.5 Update Your Account Information

You can manage user account settings in the following ways:

- Change your account information (name, title, phone, and timezone)
- Change your login password (see [Change Password](#))

To update your account information:

1. Click the **User Account** icon on the top right of the AutoEncoder status bar.



2. Select **My Account** from the drop-down menu.

The Account Information page displays.

Account Information

Prefix *(Optional)*

First Name

Lorita

Middle Name *(Optional)*

Anna

Last Name

Violeta

Suffix *(Optional)*

Title *(Optional)*

Physician

Email

lvioleta@omnicomm.com

Phone Number
(Optional)

Timezone

(UTC-05:00) Eastern Time (US & Canada) ▼

Update Account


Cancel


3. Enter/select the applicable information.
4. Click **Update Account** button.

7.6 View System Versions


You can view the database and software versions of AutoEncoder in the **About** menu.

To view system version:

1. Click the  **Main Menu** icon (top left).
2. Select **About** from the main menu.



The About AutoEncoder page displays.



This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

7.7 Return to Home Page

You can return to the home page (Dashboard) from any AutoEncoder page.

To return to the home page:

1. Locate the AutoEncoder logo on the Status Bar.




2. Click the **AutoEncoder** logo.

7.8 Remove Message

When you perform an action in AutoEncoder, a message displays until you remove it. Messages may include a confirmation (displays as green text) or warning/error (displays as red text).

To remove a message:

1. Locate the message at the top of the page.



2. Click the **X Close** icon on the right-hand side of the message.


7.9 Change Role

If assigned multiple roles, you can change roles without logging out.

Note: The System Admin user has no other roles than System Admin.

To change a role:

1. Click the role link on the top right of the AutoEncoder status bar.



2. Select the role from the drop-down menu. The new role displays in the status bar and the current page displays the associated data. (Clicking away from the drop-down retains the current role.)


Note: If your new role does not have access to the data on the page, no rows/data are displayed. The only exception is for the Coding page, the user will be returned to the Dashboard.

Table Functions


7.10 Use Table Functions


This section refers to any tabular display of data (tables/grids), such as when running a report.


On the top right-hand side of tables in AutoEncoder, you can click the following icons to perform specific actions:


 Refresh icon – updates the data in the table

 Reset Grid to Default Layout icon – resets the table to its default (original) layout

 CSV Export icon – exports the data in the table to an Excel file

 Show/Hide Columns icon – displays check boxes where to select (show) or de-select (hide) table columns

The  Information icon is located on the top left-hand side of tables in AutoEncoder. This icon provides a brief description of the table.

You can also sort and filter tables, such sorting a column in ascending/descending order, moving a column, or filtering a column on criteria. (Use the  **Reset Grid to Default Layout** icon if you want to reset the table.)


Refer to the following sections for instructions:

- [Refresh Data](#)
- [Reset Table to Default Layout](#)

- [Export Data to an Excel file](#)
- [Show Hide a Column](#)
- [Use information Icon](#)
- [Sort and Filter Tables](#)


7.11 Sort and Filter Tables

This section refers to any tabular display of data (tables), such as when running a report.

 **Important:** After sorting or filtering a table, the modification persists until you refresh the current page, navigate to another page, or log off.

You can sort and filter tables in the following ways:

- [Sort a column](#)
- [Move a column](#)
- [Show/hide a column](#)
- [Resize a column](#)
- [Filter a column on criteria](#)

 **Note:** You can [reset the table to default layout](#) if you want to restore the table to its original settings.


7.11.1 Sort a Column

You can sort columns in ascending or descending order, (except for the Actions column). The default order is ascending.

To sort a column:

1. Click a column header of a table.

The column menu displays the ▲ icon to signify it is in ascending order.




2. To reverse the sort order, click the column again. The ▼ Descending Order icon displays in the column header.
3. To revert back to the default sort order, click the column again and it displays the ▲ Ascending Order icon.



7.11.2 Move a Column


You can move a column to the desired location.

To move a column:

Drag the column header to the desired location.



The  Move icon displays when dragging the column. (A  icon displays in the background.)

Note: If you accidentally move the column outside of the table, click the  **Display/Hide** icon and select the check box of the column you moved. (See [Display Hide a Column](#) for details.)


7.11.3 Show/Hide a Column

Certain pages in AutoEncoder allow you to show/hide table columns.

To display/hide table columns:


1. From a table, such as a report, click the  **Show/Hide Columns** icon on the top right-hand side of the page.

A drop-down list displays check boxes that correspond to each table column.




2. Select/de-select the check boxes of the columns you want to display/not display. In the example above, the Description and Active columns will not display in the page and the remaining, checked columns will display in the table.


7.11.4 Filter a Column on Criteria

Certain columns contain the  **Filter** icon. When entering the text criteria, you to filter on the following:

- Equals/Not Equals
- Starts/Ends With
- Contains/Not Contains

To filter on criteria:

1. From a table, click the  **Filter** icon of the corresponding column.




2. Select the criteria from the drop-down list.
3. Enter the criteria text. In the example above, the criteria is to filter on organizations that start with "org".


 **Note:** Click the **Clear Filter** button to reset the page with no filters.

7.11.5 Resize a Column

You can resize a column to modify the space between cells.

To resize a column:

1. Using your mouse, place the cursor on right side a column boundary until the  **Resize** icon displays.





2. Drag the column boundary to the left or right until the desired width displays.


7.12 Export to Excel

You can view and/or save a report as an Excel file.




To export a report Excel:







1. Click the  **Main Menu** icon (top left).
2. Access the applicable option from the drop-down list. (After accessing System Reports, click the applicable report link.)



3. Click the  **CSV Export** icon on the right-hand side of the table.

AutoEncoder Dashboard

Select a module:   

▼ Trial	Module	Dictionary	Code	Review	Approve	Organization	Trial Group	Actions	Job Status
Trial B	 MedDRA	MedDRA Med...	793	4	3	Healthy Solutions	Trial Group AB	Reports	100% Completed
Trial B	 WHO DRG	WhoDrug B3 ...	1010	3	1	Healthy Solutions	Trial Group AB	Reports	100% Completed
Trial B	 ST2	WhoDrug B3 ...	1	4	1	Healthy Solutions	Trial Group AB	Reports	100% Completed
Trial A	 MEDRA	MedDRA Med...	0	0	0	Organization A	Trial Group AB	Reports	
Demo Study-T...	 WHO DRG	WhoDrug B3 ...	0	0	0	Healthy Solutions	Allergy	Reports	
Demo Study-T...	 MEDRA	MedDRA Med...	0	0	0	Healthy Solutions	Allergy	Reports	


13 items


4. Click the applicable button (Open or Save) to open or save the Excel file.

Dashboard (1) - Microsoft Excel

	A	B	C	D	E	F	G	H	I	J	K
1	OmniComm AutoEncoder Report: Dashboard										
2	Trial	Module	Dictionary	Code	Review	Approve	Organizati	Trial Grou	Actions	Job Status	
3	Demo Stu	MEDRA	MedDRA Mec	0	0	0	Healthy S	[object Object]			
4	Demo Stu	WHODRG	WhoDrug B3	0	0	0	Healthy S	[object Object]			
5	Demo Stu	WHODRG	WhoDrug B3	0	0	0	Anil Orgar	[object Object]			
6	Demo Stu	WHODRG	WhoDrug B3	0	0	0	OmniCom	[object Object]			
7	Demo Stu	WHODRG	WhoDrug B3	0	0	0	Healthy S	[object Object]			
8	Demo Stu	MEDRA	MedDRA Mec	0	0	0	Healthy S	[object Object]			
9	Trial A	MEDRA	MedDRA Mec	0	0	0	Organizati	[object Object]			
10	Trial B	ST2	WhoDrug B3	1	4	1	Healthy S	[object Object]		100% Completed	
11	Trial B	MedDRA	MedDRA Mec	793	4	3	Healthy S	[object Object]		100% Completed	
12	Trial B	WHODRG	WhoDrug B3	1010	3	1	Healthy S	[object Object]		100% Completed	

7.13 Use Information Icon


Table headers contain the  Information icon that provides a brief description of the feature. Hover over this icon to view the description.





7.14 Refresh Data

You can refresh a table so it displays the most up-to-date data.

To refresh data:

1. Click the  **Main Menu** icon (top left).
2. Access the applicable option from the drop-down list. (After accessing System Reports, click the applicable report link.)



3. Click the  **Refresh Data** icon on the right-hand side of the table. A data reload for the current grid occurs. The Refresh Data icon shows the progress of the refresh by turning blue and spinning in place.

AutoEncoder Dashboard

Select a module:

▲ Trial	Module	Dicti...	Code	Review	Approve	Organizat...	Trial Group	Actions	Job Status
Demo Study 1	WHODRG	WhoDrug...	2	0	0	OmniComm		Reports	100% Completed
Demo Study 1	MEDRA	MedDRA ...	20	1	0	OmniComm		Reports	
Demo Study 2	WHODRG	WhoDrug...	0	0	0	OmniComm		Reports	
Demo Study 2	MEDRA	MedDRA ...	9	0	0	OmniComm		Reports	
Demo Study Edi...	MEDRA	MedDRA ...	0	-	-	OmniComm	FLA	Reports	
P500-0717	WHODRG	WhoDrug...	0	-	-	FLA	TG A	Reports	
SDTM_DEMO	MH	MedDRA ...	9	1	1	OmniComm	FLA	Reports	
SDTM_DEMO	CPCODE	MedDRA ...	1	0	0	OmniComm	FLA	Reports	
SDTM_DEMO	MedDRA	MedDRA ...	9	1	0	OmniComm	FLA	Reports	
SDTM_DEMO	WHODRG	WhoDrug...	6	0	4	OmniComm	FLA	Reports	100% Completed
TrialA	AECODE	MedDRA ...	98	0	0	Organization A	TG A	Reports	

14 items


7.15 Reset Table to Default Layout


You can reset a table to its default (original) layout. The layout is saved automatically upon any change to the column display, such as resizing, moving, or showing/hiding a column. The current table displays these changes (adjusted for screen resolution) for the current user on reload, navigation, role change, or restart.


The Coding, Transmission Queue, Manage Synonym Terms, Reports, and Jobs pages are specific to their objects (trials, trial groups, and organizations), so changing the columns of Trial 42123_ZKS/Module AE for example, only applies to that specific page and is not reflected in other Trial/Module coding pages.

Tables settings are stored per user for the Management type pages, and also per user per module for the Coding, Transmission, reports type pages. For example, one module's Coding page can have certain settings while another module can have totally different settings.





To reset a table to the default layout:

1. Click the  **Main Menu** icon (top left).
2. Access the applicable option from the drop-down list. (After accessing System Reports, click the applicable report link.)




- Click the  **Reset Grid to Default Layout** icon on the right-hand side of the table. The current grid is reset to display the defaults, overwriting any previous changes. The new column layout, including column visibility settings are also applied to the CSV Export.

AutoEncoder Dashboard

Select a module:    

▲ Trial	Module	Dicti...	Code	Review	Approve	Organizat...	Trial Group	Actions	Job Status
Demo Study 1	WHODRG	WhoDrug...	2	0	0	OmniComm		Reports	100% Completed
Demo Study 1	MEDRA	MedDRA ...	20	1	0	OmniComm		Reports	
Demo Study 2	WHODRG	WhoDrug...	0	0	0	OmniComm		Reports	
Demo Study 2	MEDRA	MedDRA ...	9	0	0	OmniComm		Reports	
Demo Study Edi...	MEDRA	MedDRA ...	0	-	-	OmniComm	FLA	Reports	
P500-0717	WHODRG	WhoDrug...	0	-	-	FLA	TG A	Reports	
SDTM_DEMO	MH	MedDRA ...	9	1	1	OmniComm	FLA	Reports	
SDTM_DEMO	CPCODE	MedDRA ...	1	0	0	OmniComm	FLA	Reports	
SDTM_DEMO	MedDRA	MedDRA ...	9	1	0	OmniComm	FLA	Reports	
SDTM_DEMO	WHODRG	WhoDrug...	6	0	4	OmniComm	FLA	Reports	100% Completed
TrialA	AECODE	MedDRA ...	98	0	0	Organization A	TG A	Reports	

14 items


 **Note:** Filters and sorting are not reset because they are not saved as part of layout.

OBJECTS

8.0 Manage Objects

You can manage objects (trials, trial groups, organizations) in the following ways:

- View an object (see [View Trials](#), [View Trials within a Trial Group](#), [View Trial Groups](#), or [View Organizations](#))
- Add an object (see [Add a New Trial](#), [Add Trial Group](#), or [Add a New Organization](#))
- Edit an object
- Delete an object
- Edit object synonym lists
- View object reports
- Edit object settings
- Manage Mapping

 **Note:** The options are available based on user/role privileges.

8.1 Edit an Object


You can edit information of the following objects:

- Trial – name, description, associated organization/trial group, and active status
- Trial Group/Organization – name, description, and active status


To edit an object:

1. Access the applicable page. (See [View Trials](#), [View Trial Groups](#), or [View Organizations](#) for guidance.)
2. Click the object's name (in blue text) under the Name column.

Trials ⓘ New Trial			
▲ Name	Description	Active	App Connect...
Demo Study	Demonstration Study ...	✓	OmniCommMaster...
Demo Study_1	Demonstration Study ...	✓	OmniCommMaster...

 **Note:** You can also click the **Edit** link under the Actions column to edit an object.

- Enter/select the applicable information in the trial/trial group/organization Manage Trials dialog.



Manage Trials: Demo Study

Name:

Description:

Organization:

Trial Group:

☒ **Active**

Submit **Cancel**

- If you want to de-activate the object, de-select the **Active** check box. If deactivated, all objects (trials, trial groups, and organizations) owned by the trial do not automatically display on the Dashboard or Coding pages.
- Click the **Submit** button.


8.2 Delete an Object

You can delete an object (trial, trial group, organization) if you have permission.

 **Important:** Only System Users/System Admins can add or remove a trial from a trial group.

To delete an object:

- Access the applicable page. (See [View Trials](#), [View Trial Groups](#), or [View Organizations](#) for guidance.)
- Click the **Delete** link under the Action column.

Trials  New Trial						
▲ Name	Description	Active	Organization	Trial Group	Actions	
Demo Study	Demonstration Study : TG_3	✓	Healthy Solutions	Allergy	Edit	Delete Users Modules
Demo Study edit	Demonstration Study : DM...	✓	Organization A		Edit	Delete Users Modules
Demo Study-TG4	Demonstration Study : TG4	✓	OmniComm		Edit	Delete Users Modules

- Click the **Yes** button in the message. A confirmation message displays at the top of the page and the object no longer displays.

8.3 Edit Object Synonym Lists

You can edit (change inherit settings, assign, or remove) a synonym list of an object (trial, trial group, organization).

To edit object synonym lists:

1. Access the applicable page. (See [View Trials](#), [View Trial Groups](#), or [View Organizations](#) for guidance.)
2. Click the **Synonym Lists** link under the Action column.

Trials 		New Trial				
 Name	Description	Active	Organization	Trial Group	Actions    	
Trial A	Diabetes study	✓	Customer A	Trials A and B	Edit	Delete
Trial B	Allergic reactions	✓	Customer A	Trials A and B	Users	Module
Trial C	Preventing cardiac arrest	✓	Org706		Synonym Lists	Reports
					Settings	

The Synonym Lists dialog displays the dictionary (tabs), inheritance (System, Trial Group, or Organization), dictionary range of each synonym list. If the synonym list is inherited from a different object, it is indicated in the Inherited column.

Manage Trials: Synonym Lists for Trial A

Synonym Lists ⓘ

MedDRA WhoDrug B2

Assign MedDRA Synonym List ▾

Order		Synonym List	Inherited?	Disable Inherit	Dictionary Version Range	Actions
1.		Synonym List A	Yes (Organization)	<input type="checkbox"/>	MedDRA MedDRA 16.0 - MedDRA 16.1	-
2.	⬆	AEC QA2	Yes (Organization)	<input type="checkbox"/>	MedDRA MedDRA 16.0 - MedDRA 16.2	-
3.	⬆	Anil Kumar	No	-	MedDRA MedDRA 16.0 - MedDRA 16.0	Remove
4.	⬆	Test Syn List	No	-	MedDRA MedDRA 16.0 - MedDRA 16.1	Remove

Save Assignments
Cancel

Note: See [Edit Object Settings](#) for more information about the inheritance hierarchy.

3. Click the tab of the applicable dictionary.
4. If you want to disable the inheritance, click the **Disinherit** check box of the applicable synonym list. This turns off the cascade effect of inheriting a synonym list.
5. Do one of the following:
 - To assign a synonym list, select the synonym list from the **Assign WhoDrug/MedDRA Synonym List** drop-down. The list you selected is added to the bottom row of the table. To move the list to the desired row use the **Move Up** and **Move Down** buttons, then click the **Save Assignments** button.
 - To remove a synonym list, click the **Remove** link of the applicable list under the Actions column. The list is automatically removed.

8.4 View Object Reports

You can view the following reports according to each object type:

Note: You can only view reports for objects you have access to.

- Trial – Modules, Settings, Synonym Lists, Users With Roles
- Trial Group – Trials, Users With Roles
- Organization – Trials, Users, Users With Roles

To view an object report:

1. Access the applicable page. (See [View Trials](#), [View Trial Groups](#), or [View Organizations](#) for guidance.)
2. Click the **Reports** link of the applicable object under the Action column (far right).

Trials New Trial						
Name	Description	Active	Organization	Trial Group	Actions	
Trial A	Diabetes study	✓	Customer A	Trials A and B	Edit Delete Users Modules Synonym Lists	Reports Settings
Trial B	Allergic reactions	✓	Customer A	Trials A and B	Edit Delete Users Modules Synonym Lists	Reports Settings
Trial C	Preventing cardiac arrest	✓	Org706		Edit Delete Users Modules Synonym Lists	Reports Settings

The Reports page provides a link and brief description of each report.




Manage Trials: Reports for Trial A	
Trial - Modules Report	
<ul style="list-style-type: none"> • Trial - Settings Report • Trial - Synonym Lists Report • Trial - Users With Roles Report 	

3. Click the applicable report link.

The following example shows the Trial Modules report.

Trial - Modules Report

Back

Trial - Modules Report										
Module N...	Dict Name	Version	Total Terms	Coded Terms	Reviewed Terms	Approved Terms	Sent Terms			
AE	MedDRA	16.1	13	1	2	8	0			
CONMED	WhoDrug B2	Sep 2013	9	1	0	1	0			

4. Click the **Back** button to return to the previous page.

Note: You can sort and filter the report or export it to Excel. Various table options are also available. (See [Use Table Functions](#).)

8.5 Edit Object Settings

You can change the following settings of an object (trial, trial group, organization):

- MedDRA Primary Path – use primary the path only when coding against MedDRA (Yes or No)
- System Synonym Search Rules – search for synonyms by history/synonyms only or in sequence
- HistoryScope – search for the history scope stored at the trial, group, or organization level
- AutoReviewApprove – automatically review and approve terms once confirmed (True or False)

The default settings may be changed at any level (Trial, Trial Group, Organization, and System). Settings are inherited from a higher level or may be explicitly changed at a lower level. The lowest level settings take precedence. The order of precedence is shown below.

Order of Precedence

Scenario 1 (trial DOES NOT belong to a trial group)



Scenario 2 (trial DOES belong to a trial group)



This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

To edit object settings:

1. Access the applicable page. (See [View Trials](#), [View Trial Groups](#), or [View Organizations](#) for guidance.)
2. Click the **Settings** link under the Action column.

Trials 1 New Trial						
▲ Name	Description	Active	Organization	Trial Group	Actions	
Trial A	Diabetes study	✓	Customer A	Trials A and B	Edit Delete Users Modules Synonym Lists Report Settings	
Trial B	Allergic reactions	✓	Customer A	Trials A and B	Edit Delete Users Modules Synonym Lists Reports Settings	
Trial C	Preventing cardiac arrest	✓	Org706		Edit Delete Users Modules Synonym Lists Reports Settings	

The settings pop-up for the object displays. Each field is defaulted to the inherited system setting.

Manage Trials: Settings for Trial A

MedDRA Primary Path: Inherit on (Default) No
Use Primary Path only when coding against MedDRA

Synonym Search Rules: Inherit on (Trial Group) Synonyms only
Rules on searching synonyms and history matches

History Scope: Inherit on (System) History stored at Trial level
Trial history scope level

Auto Review Approve: Inherit off False
Automatically Review and Approve Terms once Confirmed

Save Settings Cancel

3. Click to the left of the **Inherit on** toggle to override the inherited value of the applicable setting(s). The toggle changes to Inherit off and the drop-down list for the setting becomes enabled.
4. Select the applicable new value(s) from the drop-down list(s).
5. Click the **Save Settings** button. A confirmation message displays, indicating that the settings were updated.

Note: If you turn **Inherit** back on, the text area shows the cascaded inherited value and its origin in the format “(Origin) Value”.

Trials

8.6 Manage Trials


You can manage trials in the following ways:


- [View trials](#)
- [Add a new trial](#)
- [Edit a trial](#)
- [Delete a trial](#)
- [Assign user to a trial](#)
- [Manage mapping](#)
- [Edit trial synonym lists](#)
- [View trial reports](#)
- [Edit trial settings](#)

8.6.1 View Trials

You can view details about each trial, including the description, status (active/inactive), organization, and trial group.

To view trials:

1. Click the  **Main Menu** icon (top left).
2. Select **Trials** from the drop-down list.



The Trials page displays the name, description, active status (check marked if active), organization, and trial group, and allowed user/role actions for each trial. (The tasks displayed in the Trials page vary depending upon permissions assigned to the active user and role).

Trials ⓘ New Trial						
▲ Name	Description	Active	Organization	Trial Group	Actions	
Demo Study	Demonstration Study : TG_3	✓	Healthy Solutions	Allergy	Edit Delete Users Modules	
Demo Study edit	Demonstration Study : DM...	✓	Organization A		Edit Delete Users Modules	
Demo Study-TG4	Demonstration Study : TG4	✓	OmniComm		Edit Delete Users Modules	


Note: You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))

8.6.2 Add a New Trial


When adding a new trial, you must specify the organization. However, the trial group is not required.

To add a new trial:

1. Access the Trials page. (See [View Trials](#) for guidance.)
2. Click the **New Trial** button.




The New Trial dialog displays.



3. Select the **Server and Database** and **Trial** from the drop-down lists.
4. Click the **Add Trial** button.

New fields display in the New Trial dialog.



Manage Trials: New Trial [X]

Name:
Demo Study

Description:
Demonstration Study : TG4

Organization:
OmniComm ▼

Trial Group:
Trial Group AB ▼

☒ Active


Submit Cancel

5. Optionally, change the **Name** and **Description**. (The Trial Name length must be at least 3 characters and no longer than 50.)
6. Select the **Organization** and **Trial Group** (not required) from the drop-down lists. (The Organization field is Read-Only unless you are a System User/System Admin.)
7. Click the **Submit** button. A confirmation message displays at the top of the page and the new organization displays in the table.


8.7 Mapping Overview

As shown in AE-TM integration process diagram, the first step to mapping is to define the TrialMaster (TM) module verbatim terms in TrialBuilder. As the user enters/updates the verbatim terms in TM, AutoEncoder (AE) matches/maps the TM module terms to AE dictionary terms.

AE-TM Integration Process



The Mapping page allows you to create mapping between TM modules and AE dictionary versions.



The screenshot shows the Mapping page interface. It is divided into three main sections: TM Module Terms (left), Match Items Area (center), and AE Dictionary Terms (right). The TM Module Terms section lists various items like WHDRG Items and CMCLSD1. The Match Items Area shows a list of items with arrows indicating matches or references to the AE Dictionary Terms. The AE Dictionary Terms section lists items like Drug Name, Generic Name, and ATC Level 1 Code. An Automap button is located in the top right corner. The Mapping Items Area is highlighted at the bottom.

- You can automatically match/map TM/AE terms via the Automap button (upper right).
- You can manually match/map TM/AE terms by dragging the terms to the corresponding column in the Match/Mapping Items area.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

Match Items Area

- Reference terms from the module are pre-defined and locked, but the verbatim term is editable. The module's verbatim term **must** have a dictionary match to be a usable module map.
- Reference terms can be either mapped/not mapped against the dictionary: Any module reference terms **not** mapped to dictionary terms display for reference purposes only. Any verbatim and module reference terms mapped to dictionary terms are usable in searches or autocoding.
- Reference fields only display in the Mapping dialog if they exist in the calling application's module definition. (If there are only 3 references then references for 4 and 5 do not display.)

Mapping Items Area


- You can map module terms only once.
- You can use dictionary terms for any number of matches.

After matching/mapping terms, AE defines what terms are written into TM and what fields are automatically populated when users enter certain terms in iCRFs listed under adverse events, concomitant medications, and medical history.

In the following example, when the TM user enters the Verbatim Term "busonal" in the CM-busonal iCRF, the Trade Name, Generic Name, Level 1 Term/Code, and Level 2 Term/Code are automatically populated with the AE dictionary terms in the Verbatim Term section of the iCRF. The matching/mapping of these fields are shown in the Mapping page example above.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

Note: TM module terms are defined in TrialBuilder in the applicable dictionary. The following figure shows how the Verbatim Term and Trade Name can be defined based on the example above. All the terms contained in the CM-busonal iCRF are listed in the Navigation Tree in the left pane. The WhoDrug Term field sets the dictionary term type for AutoEncoder's AutoMapping feature.



The screenshot displays the OmniComm interface with a Navigation Tree on the left and two configuration windows on the right.

Navigation Tree (Left Pane):

- Patients
 - New Subject
 - PAT
 - SCREEN
 - DAY 0
 - DAY 14
 - DAY 28
 - DAY 56
 - DAY 84
 - DAY 112
 - DAY 126
 - DAY 140
 - DAY 168
 - DAY 196
 - DAY 224
 - ETERM
 - ULCER ASSESSMENT
 - DRUG DISCONTINUE
 - ADVERSE EVENTS
 - CON MEDS
 - VISIT
 - Visit
 - CM
 - CRFHeader
 - CM
 - WHODRG** (highlighted with a red box)
 - CMDRVERB
 - CMDRTRNM
 - CMDECOD
 - CMCLS1
 - CMCLSCD1
 - CMCLS2
 - CMCLSCD2
 - CMCLS3
 - CMCLSCD3
 - CMCLAS
 - CMCLASCD
 - CMDRRT
 - CMDRINDC
 - CMDRCTRY

CMDRVERB - Item Configuration Window:

- Control Type: TextBox w/WhoDrug
- Data
 - Stop Word Group
 - Required: No
 - Required Text: This item is required
 - Default
 - WhoDrug Term: VERBATIM TERM** (highlighted with a red box)
- Design
 - Caption: CMDRVERB
 - OID: CMDRVERB
 - Question Text: Verbatim Term

CMDRTRNM - Item Configuration Window:

- Control Type: TextBox w/WhoDrug
- Data
 - Stop Word Group
 - Required: No
 - Required Text: This item is required
 - Default
 - WhoDrug Term: TRADE NAME
- Design
 - Caption: CMDRTRNM
 - OID: CMDRTRNM
 - Question Text: Trade Name
 - Comment
 - Question Type
 - Category: Sponsor Data
- Layout
 - Disabled: No
 - Hidden: No
 - Read Only: Yes
 - Width: 40
- Read Only

Default state of the Question upon creation. Edit must be written to overwrite default state.
- Buttons: New Item, OK, Cancel, Apply, Help

Note: See [Manage Mapping](#) to get started with the mapping process.

8.7.1 Manage Mapping

You can manage mapping in the following ways:

- View/edit a mapping between modules and a dictionary terms
- Create a mapping between modules and dictionary terms
- Upgrade dictionary version for a module's mapping
- Delete a module mapping between modules and dictionary terms

8.7.2 View/Edit a Mapping Between Modules and Dictionary Terms

You can view/edit a mapping between TrialMaster (TM) modules and AutoEncoder (AE) dictionary terms.

 **Important:** Editing a module mapping uncodes all terms and overwrites the previous mapping.

To view/edit a mapping between modules and dictionary terms:


1. Access the Trials page. (See [View Trials](#) for guidance.)
2. Click the **Modules** link of the applicable trial under the Actions column.



| Name | Description | Active | Organization | Trial Group | Actions |
|--------------|-------------------------|--------|----------------|-------------|--|
| Demo Study 1 | Demonstration Study ... | ✓ | OmniComm | | Edit Delete Users Modules Synonym Lists Reports Settings |
| Demo Study 2 | Demonstration Study ... | ✓ | OmniComm | | Edit Delete Users Modules Synonym Lists Reports Settings |
| Demo Study | Demonstration Study ... | ✓ | OmniComm | FLA | Edit Delete Users Modules Synonym Lists Reports Settings |
| P500-0717 | P500-0717 : 50024_F... | ✓ | FLA | FLA | Edit Delete Users Modules Synonym Lists Reports Settings |
| SDTM_DEMO | SDTM_DEMO : TRN_5... | ✓ | OmniComm | FLA | Edit Delete User Modules Synonym Lists Reports Settings |
| TrialA | UNIVERSE : SAE_500... | ✓ | Organization A | TG A | Edit Delete Users Modules Synonym Lists Reports Settings |

The Modules page displays each module, the dictionary/upgrade version, and mapping status (A module that was already mapped is indicated by a check mark in the Mapped column and a View link in the Mapping column).


3. Click the **View** link under the Actions column.



| Module | Mapped? | Dictionary Version | Upgrade Version | Actions |
|--------|---------|---------------------------|-----------------|--------------------|
| WHODRG | ✓ | WhoDrug C WhoDrugCMar2017 | Upgrade ... | View Delete |
| MEDRA | ✓ | WhoDrugB2Sep2013 | Upgrade ... | View Delete |

Close

The Mapping page displays the matched/mapped TM Module terms in blue in the left column of the Match/Mapping Items area. The matched/mapped AE dictionary terms display in green in the right column of the Match/Mapping Items area.



The screenshot shows the 'Mapping' interface with three main columns:

- WHODRG Items:** A list of medical terms including CMDRVERB, CMDRTRNM, CMDECOD, CMCLS1, CMCLSCD1, CMCLS2, CMCLSCD2, CMCLS3, CMCLSCD3, and CMCLAS.
- Match Items:** A central area with arrows indicating relationships. It shows 'Match' for CMDRVERB, 'Reference' for CMDRTE, CMDRINDC, and CMDRCTRY, and 'Map' for CMDRTRNM, CMDECOD, and CMCLS1.
- Matched/Mapped AE Dictionary Terms:** A list of terms mapped to the WHODRG items, including Drug Name, Generic Name, DRecNo, Seq1, Seq2, ATC Level 1 Code, ATC Level 1 Name, ATC Level 2 Code, ATC Level 2 Name, and ATC Level 3 Code.


Note: If you do not have edit permissions, the Mapping page is Read-Only for viewing.

4. To edit the mapping, do the following:

- To match/map TM module terms to AE dictionary terms, drag the applicable item from the far left column to an empty field in the first column of Match/Mapping Items area. Continue dragging items as needed. (Clicking on the **x** in the upper left corner of the item unmaps/unmatches the item.)
- To match/map AE dictionary terms to TM module terms, drag the dictionary item from the far right column to an empty field in the second column of the Match/Mapping items area. Continue dragging items as needed.

Note: For more information about mapping, refer to [Mapping Overview](#).

- Click the **Save Mapping** button. Click the **Yes** button in the Confirmation message displays, indicating that any existing mappings will be overwritten. [If the message indicates that coded terms will be invalidated (uncoded), you can click the **No** button.]



Confirmation

0 Coded terms will be invalidated if trial module is mapped. Are you sure that you want overwrite any existing mappings with changes?

Yes **No**


Note: Once a module map is saved the trial, it becomes available in the user's Dashboard if that user has access.

8.7.3 Create a Mapping Between Modules and Dictionary Terms

You can create a mapping between TrialMaster (TM) modules and AutoEncoder (AE) dictionary terms.

To create a mapping between modules and dictionary terms:

1. Access the Modules page. (See [View/Edit a Mapping Between Modules and Dictionary Terms](#) for guidance.)
2. Select the **Dictionary Version**.




Manage Trials: Modules for SDTM_DEMO

| Modules ⓘ | | | | |
|-----------|---------|---------------------------|-----------------------------|--|
| Module | Mapped? | Dictionary Version | Upgrade Version | Actions |
| WHODRG | ✓ | WhoDrug C WhoDrugCMar2017 | Upgrade ... | View
Delete |
| MEDRA | | WhoDrugB2Sep2013 ▼ | Upgrade ... | Create |

[Close](#)

3. Click the **Create** link under the Actions column.

The Mapping page displays. Terms to be coded from the calling application consist of one verbatim term to match against and up to 5 reference fields. The module's verbatim term plus 1-5 reference fields are pre-defined by TrialBuilder.



4. Click the **Automap** button to automatically match/map the TM module to AE dictionary terms.
5. Click the **Yes** button in the confirmation message to overwrite any existing matches/mappings with your new changes.

The TM module and AE dictionary terms are automatically mapped in the Mapping Items area.

Note: Automapping fields is not possible if the requested dictionary level being defined for the module's fields do not exist.

6. To match/map TM module terms to AE dictionary terms, drag the applicable item from the far left column to an empty field in the first column of Match/Mapping Items area. Continue dragging items as needed. (Clicking on the **x** in the upper left corner of the item unmaps/unmatches the item.)
7. To match/map AE dictionary terms to TM module terms, drag the dictionary item from the far right column to an empty field in the second column of the Match/Mapping items area. Continue dragging items as needed.

In the figure above, the TM module verbatim term (CMDRVERB) is matched to the AE dictionary Drug Name, the TM module country reference term (CMDRTRY) term is matched to the AE dictionary Country term. The remaining terms in the Mapping Items area are mapped except for the TM module Level 1 term (CMCLSCD2) because there is no corresponding term in the right column of the Mapping Items area.

Note: For more information about mapping, refer to [Mapping Overview](#).

8. Click the **Save Mapping** button. Click the **Yes** button in the warning message displays, indicating that any existing mappings will be overwritten.
Once a module map is saved the trial, it becomes available in the user's Dashboard if that user has access.


Note: Unmapped modules do not display in the Dashboard as they can't be coded against. If saving an existing module map, all terms are uncoded that were coded using the previous mapping.

8.7.4 Upgrade Dictionary Version for a Module's Mapping


When upgrading a dictionary version, an Impact Report indicates if any conflicts exist between the new and current dictionaries.


To upgrade the dictionary version for a module's mapping:

1. Access the Modules page. (See [View/Edit a Mapping Between Modules and Dictionary Terms](#) for guidance.)
2. Select the **Upgrade** link and select the new dictionary version.



The system navigates to the Impact Report, which displays any conflicts between the current and new dictionaries. (No information displays in the report page if there are no conflicts, as shown in the figure above.) These conflicts may consist of one of the two following values:



- Term obsolete – the term against which the user coded (directly or via synonym/history match) does not exist in the new dictionary version
 - Hierarchy changed – one or more terms above the preferred term are different in the new dictionary version (e.g., new HLT for MedDRA).
3. After reviewing the impacted terms, click the **Do Upgrade** button to continue with the Dictionary Upgrade, or **Cancel** to exit without making any changes. If you clicked **Do Upgrade**, the Modules dialog re-displays with the new dictionary version.
-  **Note:** The Upgrade Version link only displays if one or more later dictionary versions are available. Continuing with the upgrade uncodes all terms in the impact table and the dictionary version is upgraded.


8.7.5 Delete a Mapping Between Modules and Dictionary Terms

You can delete (unmap) a module entirely from its dictionary if your current user role has permissions.

 **Important:** Deleting a module mapping uncodes all terms and deletes the previous mapping.

To delete a module mapping between modules and dictionary terms:

1. Access the Modules page. (See [View/Edit a Mapping Between Modules and Dictionary Terms](#) for guidance.)
2. Click the **Modules** link of the applicable trial under the Actions column.



This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

The Modules page displays the module, the dictionary/upgrade version, and mapping status (A module that was already mapped is indicated by a check mark in the Mapped column).

Trial: Demo Study 1
✕

Modules ⓘ

| Module | Mapped? | Dictionary Version | Upgrade Version | Actions |
|--------|---------|--------------------------------|-----------------|--|
| WHODRG | ✓ | WhoDrug B2
WhoDrugB2Sep2016 | Upgrade ... | View
Delete |
| MEDRA | ✓ | MedDRA MedDRA 16.1 | Upgrade ... | View
Delete |

Close

3. Click the **Delete** link of the applicable module under the Actions column.

If it is safe to delete the module mapping, the Confirmation message indicates that the study has no coded entries associated with the module mapping.

Confirmation
✕

❓ WHODRG module mapping for trial Demo Study 1 has no coded entries associated with it and is safe to delete. Continue?

Yes
No

If existing items must be uncoded to delete the mapping, the Confirmation message indicates the number of uncoded terms within the message and on the **Yes** button. (At this point you may want to click **No** to cancel the deletion.)

Confirmation
✕

❓ There are 1 existing terms that will be uncoded for MEDRA module mapping for trial Demo Study 1 if you delete this mapping. Are you sure you want to continue?

Yes (Uncode 1)
No

4. Click the **Yes** button to confirm the deletion. Depending on how the module mapping, a success or error message displays at the top of the Trials page.

Trial Groups

8.8 Manage Trial Groups

You can manage trial groups in the following ways:


- [View trial groups](#)
- [View trials within a trial group](#)
- [Add a new trial group](#)
- [Edit a trial group](#)
- [Delete a trial group](#)
- [Assign user to a trial group](#)
- [Edit trial group synonym lists](#)
- [View trial group reports](#)
- [Edit trial group settings](#)

8.8.1 View Trial Groups

You can view details about each trial group, including the description and status (active/inactive).

To view trial groups:

1. Click the  **Main Menu** icon (top left).
2. Select **Trial Groups** from the drop-down list.



This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

The Trial Groups page displays the name, description, active status (check marked if active), and allowed user/role actions for each trial group

| Trial Groups New Trial Group | | | |
|---|--|--------|--|
| Name | Description | Active | Actions |
| ADHD | Attention Deficit Hyperactivity Disorder | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| Allergy | Test responses to allergy medications. | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FLA | FLA Center | ✓ | Edit Delete Users Trials Synonym Lists Reports |

Note: You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))

8.8.2 Add a New Trial Group

You can add a new trial group to a trial.

To add a new trial group:

1. Access the Trials Groups page. (See [View Trial Groups](#) for guidance.)
2. Click the **New Trial Group** button.

| Trial Groups New Trial Group | | | |
|---|--|--------|--|
| Name | Description | Active | Actions |
| ADHD | Attention Deficit Hyperactivity Disorder | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| Allergy | Test responses to allergy medications. | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FLA | FLA Center | ✓ | Edit Delete Users Trials Synonym Lists Reports |

The New Trial Group dialog displays.

Manage Trial Groups: New Trial Group

Name:
 Sinusitis

Description:
 Compare nose sprays

☒ Active

Submit Cancel


3. Enter the **Name** and **Description** of the trial group.
4. Click the **Submit** button. A confirmation message displays at the top of the page that a new trial group was successfully created and the new organization displays in the table.

8.8.3 View Trials within a Trial Group

You can view trials that belong to trial groups.


To view trials within a trial group:

1. Access the Trials Groups page. (See [View Trial Groups](#) for guidance.)
2. Click the **Trials** link under the Action column.



| Name | Description | Active | Actions |
|---------|--|--------|--|
| ADHD | Attention Deficit Hyperactivity Disorder | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| Allergy | Test responses to allergy medications. | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FLA | FLA Center | ✓ | Edit Delete Users Trials Synonym Lists Reports |

The Trials in Trial Group table list each trial and organization that belongs to the trial group.



| Trials in Trial Group | | |
|-----------------------|----------------------------|-------------------|
| Trial | Description | Organization |
| Demo Study | Demonstration Study : TG_3 | Healthy Solutions |
| Demo Study-TG4 | Demonstration Study : TG4 | FLA |

Note: A trial may appear in only one trial group. (To add a trial to a trial group, see [Add a New Trial](#).)

Organizations

8.9 Manage Organizations

- You can manage organizations in the following ways:
- [View organizations](#)
- [View Trials within an Organization](#)
- [Add a new organization](#)
- [Edit an organization](#)
- [Delete an organization](#)
- [Assign user to an organization](#)
- [Edit organization synonym lists](#)
- [View organization reports](#)
- [Edit organization settings](#)


This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.


8.9.1 View Organizations

You can view details about each organization, including the description and status (active/inactive).

Note: If you are not a System User/System Admin, the page only displays organizations for which your current role has at least Read permissions.


To view organizations:

1. Click the  **Main Menu** icon (top left).
2. Select **Organizations** from the drop-down list.



The Organizations page displays the name, description, active status (check marked if active), and allowed user/role actions for each organization.

| Organizations ⓘ New Organization | | | |
|---|--|--------|--|
| ▲ Name | Description | Active | Actions |
| FAU | Florida Atlantic University | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FAU Medical Research Fac... | Medical research facility at Florida Atlantic University | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FLA | FLA | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| Healthy Solutions | FL healthy solutions | ✓ | Edit Delete Users Trials Synonym Lists Reports |

 **Note:** You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))

8.9.2 Add a New Organization

You can add a new organization to a trial.

To add a new organization:

1. Access the Organizations page. (See [View Organizations](#) for guidance.)
2. Click the **New Organization** button.

| Organizations ⓘ New Organization | | | |
|---|--|--------|--|
| ▲ Name | Description | Active | Actions |
| FAU | Florida Atlantic University | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FAU Medical Research Fac... | Medical research facility at Florida Atlantic University | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| FLA | FLA | ✓ | Edit Delete Users Trials Synonym Lists Reports |
| Healthy Solutions | FL healthy solutions | ✓ | Edit Delete Users Trials Synonym Lists Reports |

The New Organization dialog displays.

Manage Organizations: New Organization
 ✕

Name:

Description:

☒ Active


3. Enter the **Name** and **Description** of the trial group.
4. Click the **Submit** button. A confirmation message displays at the top of the page that a new organization was successfully created and the organization displays in the table.

8.9.3 View Trials within an Organization

You can view trials that belong to organizations.


To view trials within an organization:

1. Access the Organizations page. (See [View Organizations](#) for guidance.)
2. Click the **Trials** link under the Action column.



| Name | Description | Active | Actions |
|----------|------------------|--------|---|
| FLA | FLA | ✓ | Edit Delete User Trials Synonym Lists Reports Settings |
| OmniComm | OmniComm Systems | ✓ | Edit Delete Users Trials Synonym Lists Reports Settings |

The Trials in Organization table list each trial and organization that belongs to the trial group.




| Trials in Organization | | |
|------------------------|-----------------------------|-------------|
| Trial | Description | Trial Group |
| Demo Study 2 | Demonstration Study : 50024 | |
| SDTM_DEMO | SDTM_DEMO : TRN_50024 | FLA |

USERS

9.0 Manage Users

There are two types of user in Autoencoder:

- Regular Users – can maintain objects owned by their organization. This user has no access to objects owned by other organizations. Regular users can only edit and grant roles in objects, manage synonym lists, and add new users that belong to his/her organization.
- System Users/System Admins – can perform operations on all objects across all organizations. Only these users can create trials, trial groups, and organizations. A trial group spans organizations, so only a System User/System Admin can create trial groups and grant permissions on trial groups.

 **Note:** Regular users can only become System Users/System Admins if they check the **System User** check box when editing the Regular User. If this user reverts back to a regular user, they have the roles that they had prior to becoming a System user. (See [Edit/Add a User](#) for more information.)


You can manage users in the following ways:


- [View users](#)
- [Edit/add a user](#)
- [Assign user role to an object](#)
- [Reset a user password](#)
- [View user reports](#)

9.1 View Users

You can view details about Regular User and System Users/System Admins.

To view users:

1. Click the  **Main Menu** icon (top left).
2. Select **Users** from the drop-down list.



The Users page displays the username, organization, title, name, email, and status. Under the System column, System Users /System Admins are indicated by a check mark. Regular users have no check mark. Check marks also indicate the user status (active/inactive, locked/unlocked).

| Users ⓘ New User | | | | | | | | | | |
|-------------------------------|--------------|--------------------|------------|-------------------|--------|--------|--------|-----------------------------------|--|--|
| Username | Organization | Title | Name | Email | Active | System | Locked | Actions | | |
| user224 | Customer A | Medical Technician | Jane Doe | user224@gmail.com | ✓ | | | Edit Roles Reset Password Reports | | |
| user261 | Customer A | Data Specialist | John Smith | user261@gmail.com | ✓ | | | Edit Roles Reset Password Reports | | |
| user267 | Customer A | Medical Reviewer | Jill Doe | user267@gmail.com | ✓ | | | Edit Roles Reset Password Reports | | |

Note: By default, users have read access to everything owned by their organization. (See [User Permissions](#) for more information.)

Note: You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions](#).)

9.2 Edit/Add a User

You can edit an existing user's information except for the user name.

If you have permissions, you can add a new user.

To edit/add a user:

1. Access the Users page. (See [View Users](#) for guidance.)
2. Click one of the following:
 - A user name (in blue text in the Users column) – to edit a user (You can also click **Edit** under the Actions column to edit a user.)
 - **New User** button – to add a new user

The Edit/New User page displays.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

Manage Users: New User

Account Information

Username:

Organization:

Active ☒

System User ☐

Password Lock ☒ No lock

User Information

Prefix (Optional)

First Name

Middle Name (Optional)

Last Name

Suffix (Optional)

Title (Optional)

Email

Phone Number (Optional)

Timezone (Optional)

3. Enter/select the applicable information:

- **Organization:** Select an organization from the drop-down menu. (This field is Read-Only for regular users; only the System User /System Admin can edit this field.)
- **Active** (default): Uncheck the Active check box to de-activate a user.
- **System User:** The System User is unchecked by default and can only be checked by the System User/System Admin.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

- **Password Lock:** "No Lock" with an icon displays for users that are not currently locked out. "Locked" with a link displays if the user has been locked out. To unlock the user, click the link. **Only a System User /System Admin or a Super User for people in his/her organization can unlock users.** Locks eventually time out and unlock automatically according to the "MinutesTillLockoutEnds" system setting. (See [View/Edit System Settings](#) for more information.) If a user is de-activated and attempts to log in, an error message displays an error: <User ID> has been deactivated in the system, contact your system admin.

Password Lock settings:

- The number of invalid password attempts is controlled by the "LockoutOnFailedLogins" setting; the default is 5
- The length of lockout period is controlled by the "MinutesTillLockoutEnds" setting; the default is 60 minutes. If this setting is zero then the user is locked out forever until a System User/System Admin or a Super User for people in his/her organization unlocks by unchecking the Locked Out check box.

Once the lockout period expires, the user can attempt logging in again. Prior to the lockout period ending, a System User /System Admin or a Super User for people in his/her organization must click the **Locked** link to allow the user to login again.

- **User Information:** The First and Last Name, and email are required; the remaining fields are optional.
4. Click one of the following buttons:
- **Update User** – to save your edits
 - **Create User** – to save the new user. The user receives three notification emails: A notification that a new account has been created, a notification with a temporary password, and a notification that the password has been changed.


A confirmation message displays at the top of the Users page, indicating that the user was successfully edited/created.

9.3 Reset a User's Password

If you are a System User/System Admin, you can reset a user's password (if the user gets locked out, etc.).

To reset a user's password:

1. Access the Users page. (See [View Users](#) for guidance.)
2. Click the **Reset Password** link under the Actions column of the applicable user.



| Username | Organization | Title | Name | Email | Active | System | Locked | Actions |
|----------|--------------|--------------------|------------|-------------------|--------|--------|--------|---|
| user224 | Customer A | Medical Technician | Jane Doe | user224@gmail.com | ✓ | | | Edit Role Reset Password Reports |
| user261 | Customer A | Data Specialist | John Smith | user261@gmail.com | ✓ | | | Edit Roles Reset Password Reports |
| user267 | Customer A | Medical Reviewer | Jill Doe | user267@gmail.com | ✓ | | | Edit Roles Reset Password Reports |


3. Click the **Yes** button in the Confirmation message. A confirmation message displays at the top of the Users page and password notification and temporary password is sent to the user's email address.

9.4 Assign User to an Object

You can assign a user to an object (organization, trial group, trial) if you are a System Administrator. Super Users may assign users within organizations and trials they are assigned to. However, Super Users cannot edit users within a trial group.


To assign a user to an object:

1. Access the Organization, Trial Groups, or Trials page. (See [View Organizations](#), [View Trial Groups](#), or [View Trials](#) for guidance.)
2. Click the **Users** link under the Actions column of the applicable organization, trial group, or trial.



| Name | Description | Active | Actions |
|--------------------|--------------------------------|--------|--|
| Allergy | Test responses to nose sprays | ✓ | Edit Delete Users Trials Synonym Lists Reports Settings |
| CUBE GROUP | CUBE | ✓ | Edit Delete Users Trials Synonym Lists Reports Settings |
| Trial Group AB | TG for testing | ✓ | Edit Delete Users Trials Synonym Lists Reports Settings |
| Trial Group1Edited | Trial Group2 descriptionEdited | | Edit Delete Users Trials Synonym Lists Reports Settings |

3. Select the user from the Assign User drop-down list.



Manage Trial Groups: Users for Allergy

Users ⓘ Assign User ▼

| User | Approver | Coder | Read | ReadOnly Monitor | Reviewer | Super User |
|-------------|-------------------------------------|--------------------------|-------------------------------------|--------------------------|--------------------------|--------------------------|
| test_admin | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| test_admin2 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| test_user | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Submit
Cancel


4. To assign roles, click the applicable check box(es) of the new user. In the example above, the test_admin user is assigned the Approver and Read-Only roles.
5. Click the **Submit** button. The user displays the table.

9.5 View User Reports

You can view user reports, such as User Role History or Coding Action reports.


To view a user report:

1. Access the Users page. (See [View Users](#) for guidance.)
2. Click the **Reports** link under the Actions column of the applicable user.




| Username | Organiz... | Title | Name | Email | Active | System | Locked | Actions |
|----------|------------|-------|-----------|-----------------|--------|--------|--------|--|
| User_1 | Org A | | user 1 | harishms@i... | ✓ | ✓ | | Edit Roles Reset Password Reports |
| testuser | OmniComm | | test user | vivek@iellip... | ✓ | | | Edit Roles Reset Password Reports |

The Manage User Reports page displays.



3. Click on a link to view a report. The report displays details about the user. For example, the User Role History Report shows a history of user actions on each object.



User_1 User - Role History Report [Back](#)

| Audit Date | By User | Action | Role | Role Description | On Object |
|-----------------------------|------------------------|--------|------|------------------|--------------------------------------|
| 19-MAR-2018 11:19:14 AM UTC | AutoEncoder Admin (#1) | Remove | Read | Read Coding | Organization: Org A (Organization A) |
| 13-MAR-2018 09:13:52 AM UTC | AutoEncoder Admin (#1) | Add | Read | Read Coding | Organization: OmniComm (OmniComm) |
| 15-MAR-2018 05:59:42 AM UTC | AutoEncoder Admin (#1) | Add | Read | Read Coding | Organization: OmniComm (OmniComm) |
| 19-MAR-2018 12:09:33 PM UTC | AutoEncoder Admin (#1) | Add | Read | Read Coding | TrialGroup: TG AB (Trial Group AB) |
| 15-MAR-2018 05:59:42 AM UTC | AutoEncoder Admin (#1) | Remove | Read | Read Coding | Organization: OmniComm (OmniComm) |
| 16-MAR-2018 04:21:52 AM UTC | AutoEncoder Admin (#1) | Add | Read | Read Coding | Organization: Org A (Organization A) |
| 16-MAR-2018 04:22:43 AM UTC | AutoEncoder Admin (#1) | Remove | Read | Read Coding | Organization: OmniComm (OmniComm) |
| 19-MAR-2018 12:04:58 PM UTC | AutoEncoder Admin (#1) | Remove | Read | Read Coding | Organization: Org A (Organization A) |
| 19-MAR-2018 12:04:50 PM UTC | AutoEncoder Admin (#1) | Remove | Read | Read Coding | Organization: Org C (Organization C) |

Note: You can sort and filter the table or export it to Excel. Various table options are also available. (See [Use Table Functions](#).)

9.6 User Permissions

The following user permissions are available.

| Permission | Description |
|------------------|--|
| Read Coding | View the coding for a module in a trial |
| Edit Coding | Perform coding for a module in a trial |
| Review Coding | Perform the Review for coding decisions |
| Approve Coding | Approve coding after review |
| Read Synonyms | View synonym lists and history matches |
| Edit Synonyms | Load and edit synonym lists and edit history matches |
| Read Mapping | Read the column mappings for a module |
| Edit Mapping | Set up and edit column mappings |
| Read Users Roles | View information about other users within their defined scope (global or organization) |
| Edit Users Roles | Can set up and edit users within their defined scope |
| Read Settings | View configuration settings for an assigned object trial, trial group, organization) |
| Edit Settings | Create and modify settings for an assigned object |
| Read Object | Read access to assigned objects (trial, trial group, organization) |
| Edit Object | Create and update access to assigned objects |

ROLES

10.0 Manage Roles

The following fixed set of roles and their permissions are pre-defined in the AutoEncoder.

Note: You cannot add a new role, nor to change the permissions for an existing role. However, you can assign roles to objects.

| Role | Permissions |
|----------------------------|--|
| Read | Views but not update assigned objects (e.g., coding for a module within a trial) |
| Coder | Performs coding and view objects that support coding, such as synonym lists |
| Reviewer | Performs a review of coding decisions |
| Approver | Approves reviewed coding decisions |
| Super User | <p>Performs all activities except the creation of objects (trials, trial groups, organizations). The Super User can edit any object owned by the user's organization.</p> <p>General security rules for non System Admin/System Users are:</p> <ul style="list-style-type: none"> • They must have a role with permissions for the object. • They must belong to that organization. <p>Note: Users should NEVER be able to manage objects outside of their organization.</p> |
| Read Only Monitor | Read only permissions, to include users |
| System Admin (System User) | <p>Automatically assigned when the user is a System Admin/System User. Provides all permissions on all objects across all organizations. This user cannot be granted to a Regular User. This user has no other roles other than System Admin/System User.</p> <p>Can perform operations on all objects across all organizations. Only these users can create trials, trial groups, and organizations. A trial group spans organizations, so only a System Admin/System User can create trial groups and grant permissions on trial groups.</p> |


You can manage roles in the following ways:


- [View roles](#)
- [Assign a role to an object](#)

10.1 View Roles

You can view the roles that have been pre-defined by AutoEncoder.

To view roles:


1. Click the  **Main Menu** icon (top left).
2. Select **Roles** from the drop-down list.



The Roles page displays the permissions for each role.

View Roles

| Roles ⓘ | | | | | | | | | | |
|------------------|-------------|-------------|---------------|----------------|---------------|---------------|---------------|---------------|------------------|------------------|
| Role | Read Coding | Edit Coding | Review Coding | Approve Coding | Read Synonyms | Edit Synonyms | Read Mappings | Edit Mappings | Read Users Roles | Edit Users Roles |
| Approver | ✓ | ✓ | | ✓ | ✓ | | ✓ | | | |
| Coder | ✓ | ✓ | | | ✓ | | ✓ | | | |
| Read | ✓ | | | | ✓ | | ✓ | | | |
| ReadOnly Monitor | ✓ | | | | ✓ | | ✓ | | ✓ | |
| Reviewer | ✓ | ✓ | ✓ | | ✓ | | ✓ | | | |
| Super User | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |
| System Admin | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ | ✓ |

 **Note:** You can show/hide table columns or export the table data to Excel. (See [Show/Hide a Column](#) or [Export to Excel](#).)

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.


10.2 Assign Role to an Object

Note: A user is created within an organization and is automatically assigned a Read role at the organization level or System Admin if System User was selected. If granted, the user acquires the corresponding permissions on all objects owned by that organization (trials, modules within trials).

You can assign a role to a user within an object (organization, trial group, trial) if your role has the correct permissions and is a member of the object's parent organization.

Object Levels

The role that is assigned to an object cascades according to the object level. For example, if a role is assigned to an organization or trial group, the role's permissions cascade down to the trial assigned to that organization or trial group. However, if a role is only assigned to a trial, the role has no permissions to the organization or trial group. (Note that a trial group is independent from the organization. An organization cannot have trial groups assigned; only trials can be assigned to trial groups.)



To assign a role to an object:

1. Access the Organization, Trial Groups, or Trials page. (See [View Organizations](#), [View Trial Groups](#), or [View Trials](#) for guidance.)
2. Click the **Users** link under the Actions column of the applicable organization, trial group, or trial.

| Trial Groups ⓘ | | New Trial Group | | | |
|--------------------|--------------------------------|-----------------|---------|--------|-------|
| ▲ Name | Description | Active | Actions | | |
| Allergy | Test responses to nose sprays | ✓ | Edit | Delete | Users |
| CUBE GROUP | CUBE | ✓ | Edit | Delete | Users |
| Trial Group AB | TG for testing | ✓ | Edit | Delete | Users |
| Trial Group1Edited | Trial Group2 descriptionEdited | | Edit | Delete | Users |

3. Select the role check box(es) of the applicable user(s).

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

Manage Trial Groups: Users for Allergy ✕

Users ⓘ Assign User ▾

| User | Approver | Coder | Read | ReadOnly Monitor | Reviewer | Super User |
|-------------|-------------------------------------|--------------------------|-------------------------------------|--------------------------|--------------------------|--------------------------|
| test_admin | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input checked="" type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| test_admin2 | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |
| test_user | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> | <input type="checkbox"/> |

Submit Cancel

- Click the **Submit** button. A confirmation message displays at the top of the page, indicating the assigned role was successful.

DICTIONARIES

11.0 Manage Dictionaries

You can manage dictionaries in the following ways:


- [View dictionaries](#)
- [Add a new dictionary version](#)
- [Delete a dictionary version](#)

11.1 View Dictionaries


You can view details about the current WHODrug and MedDRA dictionary versions.

To view dictionaries:


1. Click the  **Main Menu** icon (top left).
2. Select **Dictionaries** from the drop-down list.




The Manage Dictionaries page displays the current WHODrug and MedDRA dictionary versions and the allowed user actions.



- To view the dictionary load job status, click the **Jobs Listing** button.

The Dictionary Jobs page displays the details of each dictionary load, including the job id and status, dictionary and version, user, date submitted/scheduled, and load/format path. You can export this information to an Excel file by clicking the  **Export to Excel** icon.

| Dictionary Jobs 20 Items | | | | | | | | | |
|---------------------------------------|----------|----------------|------------|----------------|----------|----------------------|-------------------|---------------------|--|
| Job Id | Status | Status D... | Dictionary | Version | User | Date Submitt... | Message | Scheduled S... | |
| 20 | Finished | 28-MAR-2018... | MedDRA | MedDRA 18.1 | qa_admin | 28-MAR-2018 02:46... | Completed | 28-MAR-2018 02:4... | |
| 19 | Failed | 28-MAR-2018... | MedDRA | MedDRA 18.1 | qa_admin | 28-MAR-2018 02:34... | Error: Msg 486... | 28-MAR-2018 02:3... | |
| 18 | Finished | 28-MAR-2018... | MedDRA | MedDRA 18.1 | qa_admin | 28-MAR-2018 02:22... | Completed | 28-MAR-2018 02:2... | |
| 17 | Finished | 26-MAR-2018... | MedDRA | MedDRA 18.1 | AecAdmin | 26-MAR-2018 07:03... | Completed | 26-MAR-2018 07:0... | |
| 16 | Finished | 16-MAR-2018... | MedDRA | MedDRA 18.0.2 | qa_admin | 16-MAR-2018 07:31... | Completed | 16-MAR-2018 07:3... | |
| 15 | Failed | 16-MAR-2018... | MedDRA | MedDRA 18.0.2 | qa_admin | 16-MAR-2018 07:22... | Error: Msg 486... | 16-MAR-2018 07:2... | |
| 14 | Failed | 16-MAR-2018... | MedDRA | MedDRA 18.0.2 | qa_admin | 16-MAR-2018 07:15... | Error: Msg 486... | 16-MAR-2018 07:1... | |
| 13 | Finished | 16-MAR-2018... | MedDRA | MedDRA 18.0 | qa_admin | 16-MAR-2018 07:00... | Completed | 16-MAR-2018 07:0... | |
| 12 | Failed | 13-MAR-2018... | WhoDrug C | WHODrugC A... | AecAdmin | 13-MAR-2018 01:44... | Error: Msg 867... | 13-MAR-2018 01:4... | |
| 11 | Finished | 13-MAR-2018... | MedDRA | MedDRA 16.1 | qa_admin | 13-MAR-2018 01:14... | Completed | 13-MAR-2018 01:1... | |
| 10 | Finished | 13-MAR-2018... | WhoDrug B3 | WhoDrugB3 M... | qa_admin | 13-MAR-2018 12:22... | Completed | 13-MAR-2018 12:2... | |
| 9 | Finished | 12-MAR-2018... | WhoDrug B2 | WhoDrugB2 J... | qa_admin | 12-MAR-2018 06:23... | Completed | 12-MAR-2018 06:2... | |
| 8 | Finished | 12-MAR-2018... | WhoDrug B2 | WhoDrugB2 S... | qa_admin | 12-MAR-2018 06:22... | Completed | 12-MAR-2018 06:2... | |
| 7 | Finished | 12-MAR-2018... | WhoDrug B2 | WhoDrug Sep... | qa_admin | 12-MAR-2018 05:09... | Completed | 12-MAR-2018 05:0... | |
| 6 | Finished | 12-MAR-2018... | WhoDrug C3 | WhoDrugC3 M... | AecAdmin | 12-MAR-2018 02:29... | Completed | 12-MAR-2018 02:2... | |
| 5 | Finished | 12-MAR-2018... | WhoDrug C | WhoDrugC Ma... | AecAdmin | 12-MAR-2018 02:29... | Completed | 12-MAR-2018 02:2... | |
| 4 | Finished | 12-MAR-2018... | WhoDrug B3 | WhoDrugB3 M... | AecAdmin | 12-MAR-2018 02:28... | Completed | 12-MAR-2018 02:2... | |


 **Note:** You can sort and filter a table or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))

11.2 Add New Dictionary Version


You can add a dictionary version if you have permission.

To add a dictionary version:

1. Access the Manage Dictionaries page. (See [View Dictionaries](#) for guidance.)
2. Click the **New Dictionary Version** button.



The New Dictionary Version dialog displays.



The dialog box is titled 'Manage Dictionaries: New WhoDrug B3 Version'. It contains three input fields: 'Name' with the value 'March 2017', 'Load Directory' with the value 'E:\Data\AEC5Dictionary\Loads\WhoDrug\March2017', and 'WHODrug BCP Directory' with the value 'E:\Data\AEC5Dictionary\Loads\WhoDrug'. At the bottom, there are two buttons: 'Submit' and 'Cancel'. The 'Submit' button is highlighted with a red box.

3. Enter the dictionary name and load directories. The load names must be unique and the directory must exist on the SQL Server.
4. Click the **Submit** button to submit a job to load the dictionary from the source files in the load directory.

The new dictionary version displays in the Manage Dictionaries page.


- To view the load status of the dictionary version, click the **Jobs Listing** button in the top left corner of the Manage Dictionaries page.

11.3 Delete Dictionary Version


The System User/System Admin can delete the last dictionary version, as long as it is not mapped to anything.

To delete a dictionary version:

- Access the Manage Dictionaries page. (See [View Dictionaries](#) for guidance.)
- Click the **Delete** link under the Actions column.



- Click the **Yes** button when the Confirmation message displays.



A confirmation message displays at the top of the page and the dictionary version no longer displays.

SYNONYM LISTS

12.0 Manage Synonym Lists


You can manage synonym lists in the following ways:


- [View synonym lists](#)
- [Edit/add synonym lists](#)
- [Delete synonym lists](#)
- [Import/export synonym list](#)
- [Manage synonym terms](#)

12.1 View Synonym Lists

You can view details about each synonym list, including the job status, progress, type, and user.

To view synonym lists:

1. Click the  **Main Menu** icon (top left).
2. Select **Synonym Lists** from the drop-down list.



The Manage Synonym Lists page displays each synonym list name, its organization and dictionary version, and the allowed user actions.

Manage Synonym Lists

Jobs Listing

Synonym Lists ⓘ

New Synonym List

| ▲ Synonym L... | Organization | Dictionary Version ... | Actions |
|----------------|--------------|-----------------------------|---------------------------------|
| Synonym List A | FLA | MedDRA 16.0–16.1 | Edit Delete Terms Import Export |
| Synonym List B | OmniComm | WhoDrug B2 Sep 2013–Sep ... | Edit Delete Terms Import Export |
| Synonym List C | OmniComm | WhoDrug B2 Sep 2013–Sep ... | Edit Delete Terms Import Export |

3. To view the synonym list jobs, click the **Jobs Listing** button.

The Synonym List Jobs page displays the job id, status (Submitted, Processing, Failed, or Finished), progress (0 to 100%), type (import or export), list name (synonym list name), user, date submitted, message (Completed or error), and file path/name. You can export this information to an Excel file by clicking the **Export to Excel** icon.

| Synonym List Jobs 8 Items | | | | | | | | | | |
|--|----------|---------------|----------|--------|----------------|----------|---------------|-----------------|--------------|-----------------|
| Jo... | Status | Status... | Progr... | Type | List Na... | User | Date Su... | Message | Path | File |
| 8 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |
| 7 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |
| 6 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |
| 5 | Failed | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Error: Msg 7... | \\ftl-qa-... | SynListA.txt |
| 4 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA-inv... |
| 3 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |
| 2 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |
| 1 | Finished | 13-APR-201... | 100 | Import | Synonym Lis... | qa_admin | 13-APR-201... | Completed | \\ftl-qa-... | SynListA.csv |

Note: You can sort and filter the table. Various table options are also available. (See [Use Table Functions.](#))

12.2 Edit/Add Synonym Lists

You can edit/add a synonym list if you have permission.

To edit/add a synonym list:

1. Access the Manage Synonym Lists page. (See [View Synonym Lists](#) for guidance.)
2. Click one of the following:
 - A name under the Synonym List column – to edit a synonym list
 - **New Synonym List** button – to add a new synonym list

| Synonym Lists ⓘ New Synonym List | | | |
|---|--------------|-----------------------------|---|
| ▲ Synonym List | Organization | Dictionary Version ... | Actions |
| Synonym List A | OmniComm | WhoDrug B2 Sep 2013–Sep ... | Edit Delete |
| Synonym List B | (System) | WhoDrug B3 Mar 2017–Mar ... | Edit Delete |
| Synonym List C | FLA | WhoDrug B3 Mar 2017–Mar ... | Edit Delete |

The Edit/New Synonym List dialog displays.

Manage Synonym Lists: New Synonym List
✕

Name

Organization

▼

Dictionary

▼

From Version:

▼

To Version:

▼

Create Synonym List
Cancel

3. Enter/enter the synonym list name and select the organization, dictionary, and from and to versions from the drop-down lists.
4. Click one of the following buttons:
 - **Update Synonym List** – to save your edits
 - **Create Synonym List** – to save the new synonym list

The edited/new synonym list displays in the Manage Synonym Lists page.


12.3 Delete a Synonym List

You can delete a synonym list if you have permissions to edit synonyms.

Note: If the synonyms are being used by a coded term, they are automatically uncoded.


To delete a synonym list:

1. Access the Manage Synonym Lists page. (See [View Synonym Lists](#) for guidance.)
2. Click the **Delete** link of the applicable synonym list under the Actions column.



| Synonym List | Organization | Dictionary Version ... | Actions |
|----------------|--------------|-----------------------------|--|
| Synonym List A | OmniComm | WhoDrug B2 Sep 2013–Sep ... | Edit Delete Terms Import Export |
| Synonym List B | (System) | WhoDrug B3 Mar 2017–Mar ... | Edit Delete Terms Import Export |
| Synonym List C | FLA | WhoDrug B3 Mar 2017–Mar ... | Edit Delete Terms Import Export |

The Confirmation message displays, indicating no terms will be uncoded.



3. Click the **Yes** button to confirm the deletion.
 The synonym list no longer displays in the Manage Synonym Lists page.


12.4 Import/Export Synonym List Information

You can import/export synonym list information to a specific file location on the SQL Server.


Note: Only a System User/System Admin can export synonyms.

To import/export synonym list information:

1. Access the Manage Synonym Lists page. (See [View Synonym Lists](#) for guidance.)
2. Click the **Import/Export** link under the Actions column.



The Import/Export Synonym List dialog displays.



3. Enter the import/export file location.
4. Click the **Submit** button.

A confirmation message displays at the top of the page, indicating the export job has started.


The exported file contains the synonym term and up to 5 Ref fields mapped to the Primary Keys for the following:

- MedDRA LLT Code, PT Code, HLT Code, HLGT Code, SSOC Code
- WhoDrug B is the Drug Record Number, SEQ No1, SEQ No2, Country, plus all ATC Codes
- WhoDrug C is the Medicinal Product ID plus all ATC Codes

5. To view the status of the import/export job, click the **Jobs Listing** button.

The Synonym List Jobs page displays the new the import/export job, including the status (Submitted, Processing, Failed, or Finished) and progress (0 to 100%) of the import/export.

| Synonym List Jobs 13 items | | | | | |
|---|-----------|-----------------------------|----------|-----------------------------|---|
| Job Id | Status | Status Date | User | Date Submitted | Message |
| 18 | Submitted | 19-APR-2018 09:11:32 PM UTC | AecAdmin | 19-APR-2018 09:11:32 PM UTC | |
| 17 | Finished | 22-MAR-2018 03:50:00 PM UTC | qa_admin | 22-MAR-2018 03:41:41 PM UTC | Completed |
| 16 | Finished | 22-MAR-2018 03:40:00 PM UTC | qa_admin | 22-MAR-2018 03:39:12 PM UTC | Completed |
| 15 | Failed | 22-MAR-2018 03:00:01 PM UTC | qa_admin | 22-MAR-2018 02:52:10 PM UTC | Error: Msg 7330, Line 1 Cannot fetch a r... |
| 14 | Finished | 22-MAR-2018 02:50:01 PM UTC | qa_admin | 22-MAR-2018 02:45:00 PM UTC | Completed |
| 13 | Finished | 22-MAR-2018 02:40:00 PM UTC | qa_admin | 22-MAR-2018 02:38:36 PM UTC | Completed |
| 12 | Finished | 22-MAR-2018 02:30:01 PM UTC | qa_admin | 22-MAR-2018 02:29:41 PM UTC | Completed |
| 10 | Finished | 20-MAR-2018 06:20:01 PM UTC | qa_admin | 20-MAR-2018 05:58:14 PM UTC | Completed |
| 9 | Finished | 20-MAR-2018 06:10:00 PM UTC | qa_admin | 20-MAR-2018 05:57:02 PM UTC | Completed |
| 8 | Finished | 20-MAR-2018 06:00:01 PM UTC | qa_admin | 20-MAR-2018 05:56:43 PM UTC | Completed |
| 7 | Failed | 20-MAR-2018 04:50:00 PM UTC | qa_admin | 20-MAR-2018 04:41:13 PM UTC | Error: Msg 7330, Line 1 Cannot fetch a r... |
| 6 | Finished | 20-MAR-2018 04:30:00 PM UTC | qa_admin | 20-MAR-2018 04:26:28 PM UTC | Completed |
| 1 | Finished | 13-MAR-2018 08:50:04 PM UTC | qa_admin | 13-MAR-2018 08:44:07 PM UTC | Completed |

 **Note:** You can sort and filter the table. Various table options are also available. (See [Use Table Functions](#).)


12.5 Manage Synonym Terms

You can manage synonym terms in the following ways:

- [Add \(assign\) synonym terms to synonym lists](#)
- [Delete synonym terms](#)

12.5.1 Add Synonym Terms

You can add synonym terms and assign them to a synonym list.

 **Note:** Only one synonym is allowed per verbatim term.

To add a synonym term:


1. Access the Manage Synonym Lists page. (See [View Synonym Lists](#) for guidance.)
2. Click the **Terms** link of the applicable synonym list under the Actions column.

| Manage Synonym Lists Jobs Listing | | | |
|--|--------------|---------------------------|--|
| Synonym Lists New Synonym List | | | |
| ▲ Synonym List | Organization | Dictionary Version ... | Actions |
| WhoDrugB2_SynList1 | OmniComm | WhoDrug B2 WhoDrugB2Se... | Edit Delete Terms Import Export |
| WhoDrugB2_SynList2 | OmniComm | WhoDrug B2 WhoDrugB2Se... | Edit Delete Terms Import Export |
| WhoDrugB2_SynList3 | OmniComm | WhoDrug B2 WhoDrugB2Se... | Edit Delete Terms Import Export |
| WhoDrugB3_SynList | OmniComm | WhoDrug B3 WhoDrugB3Ma... | Edit Delete Terms Import Export |
| WhoDrugC_SynList | OmniComm | WhoDrug C WhoDrugCMar2... | Edit Delete Terms Import Export |
| WhoDrugC3_SynList | FLA | WhoDrug C3 WhoDrugC3Ma... | Edit Delete Terms Import Export |

3. Select the **Dictionary Version** from the drop-down list.

This information is proprietary to OmniComm Systems, Inc., and should be treated as confidential material in accordance with existing confidentiality agreements. Unauthorized use of these materials, or the data herein is strictly prohibited.

The Manage Synonym Terms page displays a list of synonym terms if any exist.




Manage Synonym Terms: WhoDrugB2_SynList1 [Back](#)

Synonym Terms [Dictionary Version](#) WhoDrugB2Sep2013 [New Synonym Term](#)

| ▲ Synonym | Reference 1 | Reference 2 | Reference 3 | Reference 4 | Refe Actions |
|-------------|-------------|-------------|-------------|-------------|------------------------|
| haloperidol | | | | | Delete |
| tylenol | | | | | Delete |

- Click the **New Synonym Term** button.

The New Synonym Terms page displays.



Manage Synonym Terms: New Synonym Term for Synonym List A

Synonym Term
aspirin

Currently Selected

| LLT Name | LLT Code | PT Name | PT Code | HLT Name | HLT Code | HLGT Name | HLGT Code | SOC Name | SOC Code |
|----------------|----------|-------------------|----------|---------------------------|----------|---------------------|-----------|-----------------|----------|
| Aspirin asthma | 10003533 | Analgesic asth... | 10063023 | Allergies to foods, fo... | 10001737 | Allergic conditi... | 10001708 | Immune syste... | 10021428 |

[Assign](#) [Cancel](#)

Term Search 3 items match: LLT Name like 'aspirin', HLT Name like 'allergies' [Clear](#)


| aspirin | LLT Code | PT Name | PT Code | allergies | HLT Code | HLGT Name | HLGT Code | SOC Name | SOC Code |
|--------------------------|----------|-------------------|----------|---------------------------|----------|---------------------|-----------|-----------------|----------|
| Asthma aspirin-sensitive | 10056818 | Analgesic asth... | 10063023 | Allergies to foods, fo... | 10001737 | Allergic conditi... | 10001708 | Immune syste... | 10021428 |
| Aspirin asthma | 10003533 | Analgesic asth... | 10063023 | Allergies to foods, fo... | 10001737 | Allergic conditi... | 10001708 | Immune syste... | 10021428 |
| Aspirin-sensitive asthma | 10003534 | Analgesic asth... | 10063023 | Allergies to foods, fo... | 10001737 | Allergic conditi... | 10001708 | Immune syste... | 10021428 |

- Enter the **Synonym Term** and **Reference Terms** (if applicable) in the top portion of the page. (You may enter up to 5 reference terms.)
- Type the applicable search parameters in the corresponding fields and tab out or press the ENTER key on your keyboard. For example, to search for busonal terms relating to respiratory, type "busonal" in the Drug Name field and "respiratory" in the ATC Level 1 Name field (as shown in the example above).

Note: You can clear the search by clicking on the **Clear** button.

- Choose one of the search results by clicking on the item so that it becomes highlighted.
- Click the **Assign** button in the right middle portion of the page.

The Manage Synonym Terms page displays a confirmation message and the new synonym term.



Synonym term created.

Manage Synonym Terms: WhoDrugB2_SynList1 [Back](#)

Synonym Terms [Dictionary Version](#) WhoDrugB2Sep2013 [New Synonym Term](#)

| ▲ Synonym | Reference 1 | Reference 2 | Reference 3 | Reference 4 | Refe Actions |
|-------------|-------------|-------------|-------------|-------------|------------------------|
| busanol | US | | | | Delete |
| haloperidol | | | | | Delete |
| tylenol | | | | | Delete |

Note: You can sort and filter the table. Various table options are also available. (See [Use Table Functions.](#))

12.5.2 Delete Synonym Terms

You can delete synonym terms that have been assign them to a synonym list.


To delete a synonym term:

1. Access the Manage Synonym Lists page. (See [View Synonym Lists](#) for guidance.)
2. Click the **Terms** link of the applicable synonym list under the Actions column.




The Manage Synonym Terms page displays.

3. Click the **Delete** link of the applicable synonym term under the Actions column.



The Confirmation message displays, indicating no terms will be uncoded.



4. Click the **Yes** button to confirm the deletion.

The synonym list no longer displays in the Manage Synonym Terms page.


SYSTEM REPORTS


13.0 View System Reports

You can view the following system reports if you are a System Use/System Admin:


- Autoencoder Job Queue – AutoEncoder Job Statuses
- AppLog – Last 30 Days – last 30 days of the AppLog in descending order
- Error AppLog Last 30 Days – last 30 days of Errors in the AppLog in descending order
- Trials Currently Processed by AutoEncoder Service – current trials being processed by AECWSRV connected to that Autoencoder

To view a system report:

1. Click the  **Main Menu** icon (top left).
2. Select **System Reports** from the drop-down list.



The System Reports page displays report links.



3. Click on the applicable report link.

The report displays. For example, the AutoEncoder Job Queue report (below).

AutoEncoder Job Queue [Back](#)

AutoEncoder Job Queue 20 Items

| Job Id | Status | Status Date | Progress | Trial Name | Module Na... | Description | User | Date Sub... | Message |
|--------|----------|------------------|----------|--------------|--------------|-------------------|----------|------------------|-----------|
| 20 | Finished | 16-JAN-2018 0... | 100 | Demo Study_1 | AEDetails(1) | Demonstration ... | AecAdmin | 16-JAN-2018 0... | Completed |
| 19 | Finished | 16-JAN-2018 0... | 100 | Demo Study_1 | AEDetails(1) | Demonstration ... | AecAdmin | 16-JAN-2018 0... | Completed |
| 18 | Finished | 16-JAN-2018 0... | 100 | Demo Study_1 | AEDetails(1) | Demonstration ... | AecAdmin | 16-JAN-2018 0... | Completed |
| 17 | Finished | 15-JAN-2018 0... | 100 | UNIVERSE | AECODE | UNIVERSE : SA... | AecAdmin | 15-JAN-2018 0... | Completed |
| 16 | Finished | 15-JAN-2018 0... | 100 | UNIVERSE | AECODE | UNIVERSE : SA... | AecAdmin | 15-JAN-2018 0... | Completed |


- Note:** You can sort and filter the report or export it to Excel. Various table options are also available. (See [Use Table Functions.](#))


SYSTEM SETTINGS

14.0 View/Edit System Settings

You can view/edit system-wide settings. You can only edit these settings if you are a System User/System Admin.

To view/edit system-wide settings:

1. Click the  **Main Menu** icon (top left).
2. Select **System Settings** from the drop-down list.



The System Settings page displays the default value for each setting. These settings cascade to objects (trials, trial groups, organizations). Settings that are set at the object level override any system settings. (See [Edit Object Settings](#) for more information.)

System Settings

Env Name:

Environment Name

Maint Mode:

Maintenance mode toggle

System Message:

System Message

MedDRA Primary Path:

Use Primary Path only when coding against MedDRA

Synonym Search Rules:

Rules on searching synonyms and history matches

Time Format12h:

12h DateTime format

Time Format24h:

24h DateTime format

Lockout On Failed Logins:

Number of failed logins before Locking the user out

Minutes Till Lockout Ends:

Number of minutes until Lockout resets and user can try to login again

Access Token Expire Time Span:


Access Token Expire TimeSpan in minutes

Allow Insecure Http:

Allow Insecure Http for Access Token

3. To edit a setting, enter the information or select it from the drop-down lists. (Read-only fields are disabled and greyed out.)

The description is provided below each field in the System Settings page. The following provides more detail about these settings:

- Env (Environment) Name – the name you choose for this installation of AutoEncoder. This is used as an identifying label in places like the About page (click  > **About**)
- Maint (Maintenance) mode – used for maintenance issues. If set to True, regular users are blocked from logging in. System Users (System Admins) can still login.
- System Message – text that displays in the About AutoEncoder page. This also shows on the About page, and is used in addition to the default Maintenance Mode message shown to users on the Login page.
- MedDRA Primary Path – if this is set to 'Yes', then for all MedDRA dictionaries, only entries marked as Primary Path are used for coding across the entire application.
- System Search Rules – allows the administrator to set the priority of synonyms and histories when making autocoding decisions. You can select to use History only, Synonyms only, or use both with differing priorities (History then Synonyms, Synonyms then History). The settings are as follows:
 - None – only display results from dictionary
 - History only – only display History first and then dictionary results
 - Synonyms only – only display synonyms first and then dictionary results
 - Synonyms then History – display synonyms first, then history and lastly dictionary results
 - History then Synonyms – display history first, then synonyms and lastly dictionary results
- Time Format12h/24h – date format used across the application for users who have selected to display twelve-hour/twenty-four-hour format for dates and/or times. *For example, "User A" could have the 12-hour format and "User B" could have the 24-hour format.*
- Lockout on Failed Logins – an integer ranging from 1 to 15 that indicates how many unsuccessful login attempts a user can make before the system automatically locks that user out.
- Minutes Till Lockout Ends – an integer ranging from 0 to 1440 that indicates that a user who has been locked out due to too many unsuccessful login attempts must wait this amount of time before the system unlocks their account. If set to 0, only a System User/System Admin or a Super User for people in his/her organization can unlock.
- Access Token Expire Time Span – indicates the time-to-live for the AutoEncoder API access token, and may range from 10 to 1440 minutes.
- Allow Insecure HTTP – allows insecure HTTP for access token. If set to False, then HTTPS must be used. Therefore on a production system it would be false to enforce HTTPS links to the API.
- Username Reg Ex – validation pattern used for allowable user names in the application. It is not alterable in the current version of AutoEncoder.
- Username Warning – message shown when a proposed user name does not match the Username Reg Ex pattern.
- Password Reg Ex – validation pattern used for allowable user passwords in the application. It is not alterable in the current version of AutoEncoder.
- Password Warning – message shown when a proposed user password does not match the Password Reg Ex pattern.

- History Scope – application-wide setting that determines on what level trials keep their coding history: unshared (kept at their own trial level), Organization (all trials store their history in their organization, which is then available to all trials under that organization), or Group (all trials store their history in their trial group, if they have one, which is then available to all trials under that trial group). This "global" application scope can be overridden at the Organization, Trial Group or individual Trial levels.
 - Auto Review Approve – If set to TRUE, then all trials in the application (unless overridden) only present one grid on their Coding Page, the 'Coding Grid'. Any entries in the coding grid are sent immediately to the Transmission Queue on confirmation, rather than having to pass through the Review and Approval stages.
 - Days Till Password Expires – User passwords last this many days before expiring. This is an integer ranging from 0 to 365. Setting the password expiry to 0 disables all expiration.
 - Password Days Warning: This is the number of days prior to password expiration that the user will start being warned of the approaching expiration and is encouraged to change their password. It is an integer ranging from 1 to 10.
 - Pass (Password) Gen (Generation) Length: The system generates random passwords for new users and users who have forgotten their passwords. This setting is used to determine how many characters the random password should use. It is not alterable in this version of AutoEncoder.
 - Password Generation Use Upper Alpha: This Boolean value determines whether or not the random password generator uses uppercase characters. It is not alterable in this version of AutoEncoder.
 - Password Generation Use Lower Alpha: This Boolean value determines whether or not the random password generator uses lowercase characters. It is not alterable in this version of AutoEncoder.
 - Password Generation Use Numbers: This Boolean value determines whether or not the random password generator uses numeric characters. It is not alterable in this version of AutoEncoder.
 - Password Generation Use Punctuation: This Boolean value determines whether or not the random password generator uses any available non-alphanumeric characters. It is not alterable in this version of AutoEncoder.
 - Include Time Zone Information: This Boolean setting determines whether or not specific time zone settings are included in timestamps used and displayed by the application.
4. Click the **Save Settings** button on the bottom of the System Settings page.